

# Product Datasheet

## PDCD6 Antibody - BSA Free

### NBP2-94061-0.1ml

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-94061](http://www.novusbio.com/NBP2-94061)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-94061](http://www.novusbio.com/reviews/destination/NBP2-94061)



**NBP2-94061-0.1ml**

PDCD6 Antibody - BSA Free

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at -20C. Avoid freeze-thaw cycles.                                                    |
| Clonality           | Polyclonal                                                                                  |
| Preservative        | 0.02% Sodium Azide                                                                          |
| Isotype             | IgG                                                                                         |
| Purity              | Affinity purified                                                                           |
| Buffer              | PBS (pH 7.3), 50% glycerol                                                                  |


| Product Description |                                                                                                                                                                                                                                                                                                                                 |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Knockout (KO) Validated Rabbit PDCD6 Antibody - BSA Free (NBP2-94061) is a polyclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                            |
| Host                | Rabbit                                                                                                                                                                                                                                                                                                                          |
| Gene ID             | 10016                                                                                                                                                                                                                                                                                                                           |
| Gene Symbol         | PDCD6                                                                                                                                                                                                                                                                                                                           |
| Species             | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                               |
| Immunogen           | Recombinant fusion protein containing a sequence corresponding to amino acids 1-191 of human PDCD6 (NP_037364.1).<br>MAAYSYPGPGAGPGPAAGAALPDQSFLWNVFQRVDKDRSGVISDTTELQQAL<br>SNGTWTPFNPVTVRSIISMFDRENKAGVNFSEFTGVWKYITDWQNVFRTYDRD<br>NSGMIDKNELKQALSGFGYRLSDQFHDILIRKFDRQGRGQIAFDDFIQGCIVLQR<br>LTDIFRRYDTDQDGWIVSYEQYLSMVFSIV |

| Product Application Details |                                                                                   |
|-----------------------------|-----------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunoprecipitation, Knockout Validated                             |
| Recommended Dilutions       | Western Blot 1:500 - 1:1000, Immunoprecipitation 1:50 - 1:200, Knockout Validated |




## Images

Immunoprecipitation-NBP2-94061-NBP2-94061-Analysis of 200 ug extracts of HepG2 cells, using 3 ug PDCD6 antibody. Western blot was performed from the immunoprecipitation using PDCD6 antibody at a dilution of 1:1000.



Western Blot: Rabbit Polyclonal PDCD6 Antibody [NBP2-94061] - Analysis of extracts from wild type (WT) and PDCD6 knockout (KO) 293T(KO) cells, using PDCD6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-94061-0.1ml**

---

|                   |                                                     |
|-------------------|-----------------------------------------------------|
| HAF008            | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160            | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891        | Rabbit IgG Isotype Control                          |
| NBP1-51048-0.05mg | Recombinant Human PDCD6 His Protein                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-94061](http://www.novusbio.com/reviews/submit/NBP2-94061)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

