

Product Datasheet

PARG Antibody - Azide and BSA Free **NBP2-94199-0.1ml**

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-94199

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-94199



NBP2-94199-0.1ml

PARG Antibody - Azide and BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Proclin 300
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol


Product Description	
Description	Novus Biologicals Rabbit PARG Antibody - Azide and BSA Free (NBP2-94199) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	8505
Gene Symbol	PARG
Species	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 777-976 of human PARG (NP_003622.2). HNECLIITGTEQYSEYTGYAETRWRSRSHEDGSRDDWQRRCTEIVAILDALHF RRYLDQFVPEKMRRELNKAYCGFLRPGVSSENLSAVATGNWGCAGFGGDAR LKALIQILAAAAAERDVVYFTFGDSELMRDIYSMHIFLTERKLTVGDVYKLLLRYY NEECRNCSTPGPDIKLYPFIYHAVESCAETADHSGQRTGT

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:100 - 1:500



Images

Western Blot: PARG Antibody [NBP2-94199] - HCC1937 human breast cancer cells untreated (1) and treated with PARP inhibitor (2) for 96 h. 3 PARG bands 95 kDa can be seen in both cases. Primary antibody dilution 1:1000 (overnight, 4C). Image from verified customer review.



Western Blot: PARG Antibody - Azide and BSA Free [PARG] - Western blot analysis of lysates from PC-3 cells, using PARG Rabbit pAb at 1:400 dilution.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution.

Lysates/proteins: 25ug per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit .

Exposure time: 30s.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-94199-0.1ml

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00008505-Q01-10ug	Recombinant Human PARG GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-94199

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

