

# Product Datasheet

## TLE1 Antibody [PE] NBP2-98893PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-98893PE](http://www.novusbio.com/NBP2-98893PE)

Updated 10/28/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-98893PE](http://www.novusbio.com/reviews/destination/NBP2-98893PE)



**NBP2-98893PE**

TLE1 Antibody [PE]

| Product Information         |                                                                                                           |
|-----------------------------|-----------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                    |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.               |
| Storage                     | Store at 4C in the dark.                                                                                  |
| Clonality                   | Polyclonal                                                                                                |
| Preservative                | 0.05% Sodium Azide                                                                                        |
| Isotype                     | IgG                                                                                                       |
| Conjugate                   | PE                                                                                                        |
| Purity                      | Antigen and protein A Affinity-purified                                                                   |
| Buffer                      | PBS                                                                                                       |
| Product Description         |                                                                                                           |
| Host                        | Rabbit                                                                                                    |
| Gene ID                     | 7088                                                                                                      |
| Gene Symbol                 | TLE1                                                                                                      |
| Species                     | Human                                                                                                     |
| Immunogen                   | Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human TLE1. |
| Product Application Details |                                                                                                           |
| Applications                | Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                       |
| Recommended Dilutions       | Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                       |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-98893PE**

---

|                    |                                             |
|--------------------|---------------------------------------------|
| NBP2-24983         | Rabbit IgG Isotype Control [PE]             |
| H00007088-P01-10ug | Recombinant Human TLE1 GST (N-Term) Protein |
| DY453-05           | CXCL1/GRO alpha/KC/CINC-1 [Biotin]          |
| NBL1-16944         | TLE1 Overexpression Lysate                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-98893PE](http://www.novusbio.com/reviews/submit/NBP2-98893PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

