

Product Datasheet

Human LRP-6 ELISA Kit (Colorimetric) NBP3-00411

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-00411

Updated 1/8/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-00411




NBP3-00411**Human LRP-6 ELISA Kit (Colorimetric)**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Gene ID	4040
Gene Symbol	LRP6
Species	Human
Specificity/Sensitivity	This kit recognizes Human LRP-6 in samples. No significant cross-reactivity or interference between Human LRP-6 and analogues was observed.
Kit Components	Micro ELISA Plate (Dismountable), Reference Standard, Concentrated Biotinylated Detection Ab (100x), Concentrated HRP Conjugate (100x), Sample Diluent, Biotinylated Detection Ab Diluent, HRP Conjugate Diluent, Concentrated Wash Buffer (25x), Substrate Reagent, Stop Solution, Plate Sealer, Product Manual
Standard Curve Range	0.31 - 20 ng/mL
Sensitivity	0.19 ng/mL
Inter Assay	10.05%
Intra Assay	8.73%
Assay Type	Sandwich-ELISA
Spiking Recovery	87 - 106%
Suitable Sample Type	Serum, plasma, cell supernatant, and other biological fluids
Sample Volume	100 ul
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA




Images

ELISA: Human LRP-6 ELISA Kit (Colorimetric) [NBP3-00411] - Samples were spiked with high concentrations of Human LRP-6 and diluted with Reference Standard & Sample Diluent to produce samples with values within the range of the assay.



ELISA: Human LRP-6 ELISA Kit (Colorimetric) [NBP3-00411] - Standard Curve Reference





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-00411

5036-WN-010	Wnt-3a [Unconjugated]
MAB1505	LRP-6 Antibody (255302) [Unconjugated]
AF1589	SOST/Sclerostin Antibody [Unconjugated]
5439-DK-010	Dkk-1 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-00411

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

