

Product Datasheet

GCC185 Antibody - Azide and BSA Free **NBP3-03436-100ul**

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-03436

Updated 11/7/2023 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-03436



NBP3-03436-100ul

GCC185 Antibody - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Affinity purified
Buffer	PBS with 50% glycerol, pH7.3.

Product Description


Host	Rabbit
Gene ID	9648
Gene Symbol	GCC2
Species	Human, Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GCC185 (NP_852118.1). MEDLVQDGVASPATPGTGKSKLETLPKEDLIKFAKKQMMLIQAKASRCTELEKE IEELRSKPVTEGTGDIKALTERLDALLLEKAETEQQCLSLKKENIKMKQEVEDSV TKMGDAHKELEQSHINYVKEIENLKNELMAVRSKYSEDKANLQKQLEEAMNTQ LLESEQLKFQNNSEDNVKKLQEEIEKIRPGFEEQILYLQKQLDATTDEKKETVTQ LQNIIEANSQHYQKNINSLQEELLQLKAHQEEVKELMCQIEASAKEHEAEINKLN ELKENLVKQCEASEKNIQKKYECLE

Product Application Details


Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:2000, Immunohistochemistry 1:50-1:200, Immunohistochemistry-Paraffin

Images


Western Blot: GCC185 Antibody [NBP3-03436] - Analysis of extracts of various cell lines, using GCC185 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit.




Immunohistochemistry-Paraffin: GCC185 Antibody [NBP3-03436] - Paraffin-embedded mouse spleen using GCC185 antibody at dilution of 1:150 (40x lens).




Immunohistochemistry-Paraffin: GCC185 Antibody [NBP3-03436] - Paraffin-embedded human colon using GCC185 antibody at dilution of 1:150 (40x lens).



Immunohistochemistry-Paraffin: GCC185 Antibody [NBP3-03436] - Paraffin-embedded human stomach using GCC185 antibody at dilution of 1:150 (40x lens).



Immunohistochemistry-Paraffin: GCC185 Antibody [NBP3-03436] - Paraffin-embedded human appendix using GCC185 antibody at dilution of 1:150 (40x lens).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-03436-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-89788PEP	GCC185 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-03436

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

