

# Product Datasheet

## EXOC2 Antibody - BSA Free

### NBP3-04535-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-04535](http://www.novusbio.com/NBP3-04535)

Updated 4/30/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-04535](http://www.novusbio.com/reviews/destination/NBP3-04535)



**NBP3-04535-100ul**


EXOC2 Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                  |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                                                                                                                                                                                                                                           |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                      |
| Storage                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                         |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                                       |
| Preservative                | 0.02% Sodium Azide                                                                                                                                                                                                                                                                                                                                               |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                              |
| Purity                      | Affinity purified                                                                                                                                                                                                                                                                                                                                                |
| Buffer                      | PBS (pH 7.3), 50% glycerol                                                                                                                                                                                                                                                                                                                                       |
| Target Molecular Weight     | 104 kDa                                                                                                                                                                                                                                                                                                                                                          |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                  |
| Host                        | Rabbit                                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 55770                                                                                                                                                                                                                                                                                                                                                            |
| Gene Symbol                 | EXOC2                                                                                                                                                                                                                                                                                                                                                            |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                |
| Immunogen                   | Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human EXOC2 (NP_060773.3).<br>MSRSRQPPLVTGISPNEGIPWTKVTIRGENLGTGPTDLIGLTICGHNCLLTAEW<br>MSASKIVCRVGQAKNDKGDIIVTTKSGGRGTSTVSFKLLKPEKIGILDQSAVWV<br>DEMNYDDMRTDRNKGIPPLSLRPANPLGIEIEKSKFSQKDLEMLFHGMSADFTS<br>ENFSAAWYLIENHSNTSFEQLKMAVTNLKRQANKKSEGLAYVKGGLSTFFEA<br>QDALS |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                  |
| Applications                | Western Blot, Immunoprecipitation                                                                                                                                                                                                                                                                                                                                |
| Recommended Dilutions       | Western Blot 1:500 - 1:2000, Immunoprecipitation 1:50 - 1:200                                                                                                                                                                                                                                                                                                    |




## Images

Western Blot: EXOC2 Antibody [NBP3-04535] - Analysis of extracts of various cell lines, using ExOC2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit



Immunoprecipitation: EXOC2 Antibody [NBP3-04535] - Immunoprecipitation analysis of 300ug extracts of HeLa cells using 3ug EXOC2 antibody (NBP3-04535). Western blot was performed from the immunoprecipitate using EXOC2 antibody (NBP3-04535) at a dilution of 1:1000.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-04535](http://www.novusbio.com/reviews/submit/NBP3-04535)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

