

Product Datasheet

HAUS3 Antibody - BSA Free

NBP3-05007-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-05007

Updated 11/7/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-05007



NBP3-05007-100ul

HAUS3 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS with 50% glycerol, pH7.3.

Product Description


Host	Rabbit
Gene ID	79441
Gene Symbol	HAUS3
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 384-603 of human HAUS3 (NP_078787.2). FELLQLSYEIELRKHRDIYRQLENLVQELSQSNMMLYKQLEMLTDPSVSQQINP RNTIDTKDYSTHRLYQVLEGENKKKELFLTHGNLEEVAEKQKQNISLVQDQLAV SAQEHSSFFLSKRNKDMDMLCDTLYQGGNQLLLSDQELTEQFHKVESQLNKLN HLLTDILADVKTARKTLANNKLHQMEREFYVYFLKDEDYLDKDIVENLETQSKIKA VSLED

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 1:500-1:2000

Images

Western Blot: HAUS3 Antibody [NBP3-05007] - Analysis of extracts of various cell lines, using HAUS3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit

**Publications**

Ustinova K, Ruhnaw F, Gili M, Surrey T Microtubule binding of the human augmin complex is directly controlled by importins and Ran-GTP Journal of cell science 2023-06-15 [PMID: 37357828] (WB, Insect)

Details:
1:1000 WB dilution



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NBP3-05007-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-57796PEP	HAUS3 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-05007

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

