

Product Datasheet

ZNF177 Antibody - Azide and BSA Free NBP3-05074-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-05074

Updated 11/7/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-05074



NBP3-05074-100ul

ZNF177 Antibody - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Affinity purified
Buffer	PBS with 50% glycerol, pH7.3.

Product Description


Host	Rabbit
Gene ID	7730
Gene Symbol	ZNF177
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-110 of human ZNF177 (NP_003442.2). LASVG YQLCRHSLISKVDQEQLKTDERGILQGDCADWETQLKPKDTIAMQNIPG GKTSNGI

Product Application Details


Applications	Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry 1:50-1:200, Immunocytochemistry/ Immunofluorescence 1:50-1:200, Immunohistochemistry-Paraffin

Images

Western Blot: ZNF177 Antibody [NBP3-05074] - Paraffin-embedded mouse brain using ZNF177 antibody at dilution of 1:100 (40x lens).




Immunocytochemistry/Immunofluorescence: ZNF177 Antibody [NBP3-05074] - Analysis of U-2 OS cells using ZNF177 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunohistochemistry-Paraffin: ZNF177 Antibody [NBP3-05074] - Paraffin-embedded mouse kidney using ZNF177 antibody at dilution of 1:100 (40x lens).




Immunocytochemistry/Immunofluorescence: ZNF177 Antibody [NBP3-05074] - Analysis of C6 cells using ZNF177 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunocytochemistry/Immunofluorescence: ZNF177 Antibody [NBP3-05074] - Analysis of NIH-3T3 cells using ZNF177 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunohistochemistry-Paraffin: ZNF177 Antibody [NBP3-05074] - Paraffin-embedded rat pancreas using ZNF177 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry-Paraffin: ZNF177 Antibody [NBP3-05074] - Paraffin-embedded mouse pancreas using ZNF177 antibody at dilution of 1:100 (40x lens).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-05074-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00007730-P01-10ug	Recombinant Human ZNF177 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-05074

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

