

Product Datasheet

MAN2B1 Antibody - BSA Free

NBP3-05138-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-05138

Updated 1/30/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-05138



NBP3-05138-100ul

MAN2B1 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol
Target Molecular Weight	114 kDa

Product Description


Host	Rabbit
Gene ID	4125
Gene Symbol	MAN2B1
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 640-850 of human MAN2B1 (NP_000519.2). TFFWYNASIGDNESDQASGAYIFRPNQKPLPVSRWAQIHLVKTPLVQEVHQN FSAWCSQVVRLYPGQRHLELEWSVGPIVGDWKGKEVISRFDTPLETKGRFYT DSNGREILERRRDYRPTWKLNQTEPVAGNYYPVNTRIYITDGNMQLTVLTDRS QGGSSLRDGSLELMVHRLLKDDGRGVSEPLMENGSGAWVRGRHLVLLDTA Q

Product Application Details

Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500-1:2000, ELISA

Images

Western Blot: MAN2B1 Antibody [NBP3-05138] - Analysis of extracts of mouse liver, using MAN2B1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit



Publications

Hana Yuasa, Naho Morino, Takumi Wagatsuma, Masayuki Munekane, Sachiko Ueda, Mayu Matsunaga, Yasuo Uchida, Takane Katayama, Toshihiko Katoh, Taiho Kambe ZNT5-6 and ZNT7 play an integral role in protein N - glycosylation by supplying Zn 2+ to Golgi α -mannosidase II The Journal of Biological Chemistry 2024-05-16 [PMID: 38762179]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-05138-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-81326PEP	MAN2B1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-05138

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

