

# Product Datasheet

## **SERPINB1 Antibody (04) - Azide and BSA Free** **NBP3-05846-100ul**

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-05846](http://www.novusbio.com/NBP3-05846)

Updated 2/26/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-05846](http://www.novusbio.com/reviews/destination/NBP3-05846)



**NBP3-05846-100ul**

SERPINB1 Antibody (04) - Azide and BSA Free

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 100 ul                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 04                                                                                          |
| Preservative        | No Preservative                                                                             |
| Isotype             | IgG1                                                                                        |
| Purity              | Protein A purified                                                                          |
| Buffer              | 0.2 um filtered solution in PBS                                                             |

| Product Description     |                                                                                                                                                                                                                             |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description             | This antibody can be stored at 2C-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Avoid repeated freeze-thaw cycles. |
| Host                    | Mouse                                                                                                                                                                                                                       |
| Gene ID                 | 1992                                                                                                                                                                                                                        |
| Gene Symbol             | SERPINB1                                                                                                                                                                                                                    |
| Species                 | Human                                                                                                                                                                                                                       |
| Reactivity Notes        | No cross-reactivity in ELISA with:<br/><br/>Insect cell lysate                                                                                                                                                              |
| Specificity/Sensitivity | No cross-reactivity in ELISA with:<br><br>Insect cell lysate                                                                                                                                                                |
| Immunogen               | This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human SERPINB1 (Uniprot#: P30740; Met 1-Pro 379).          |

| Product Application Details |                     |
|-----------------------------|---------------------|
| Applications                | ELISA               |
| Recommended Dilutions       | ELISA 1:1000-1:2000 |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-05846-100ul**

---

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B            | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97005-0.5mg   | Mouse IgG1 Isotype Control (MG1)                        |
| H00001992-P01-10ug | Recombinant Human SERPINB1 GST (N-Term) Protein         |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-05846](http://www.novusbio.com/reviews/submit/NBP3-05846)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

