

# Product Datasheet

## NAIF1 Antibody - BSA Free

### NBP3-06197-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-06197](http://www.novusbio.com/NBP3-06197)

Updated 2/26/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-06197](http://www.novusbio.com/reviews/destination/NBP3-06197)



**NBP3-06197-100ul**

NAIF1 Antibody - BSA Free

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ul                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | 0.03% Proclin 300                                                                           |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Purity</b>        | Antigen and protein A Affinity-purified                                                     |
| <b>Buffer</b>        | PBS                                                                                         |

**Product Description**


|                    |                                                                                                                                                                                                                             |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | This antibody can be stored at 2C-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Avoid repeated freeze-thaw cycles. |
| <b>Host</b>        | Rabbit                                                                                                                                                                                                                      |
| <b>Gene ID</b>     | 203245                                                                                                                                                                                                                      |
| <b>Gene Symbol</b> | NAIF1                                                                                                                                                                                                                       |
| <b>Species</b>     | Human                                                                                                                                                                                                                       |
| <b>Immunogen</b>   | Produced in rabbits immunized with E. coli-derived Human NAIF1 fragment.                                                                                                                                                    |

**Product Application Details**

|                              |                                                     |
|------------------------------|-----------------------------------------------------|
| <b>Applications</b>          | Immunocytochemistry/ Immunofluorescence             |
| <b>Recommended Dilutions</b> | Immunocytochemistry/ Immunofluorescence 1:100-1:500 |

**Images**

Immunocytochemistry/Immunofluorescence: NAIF1 Antibody [NBP3-06197] - Immunofluorescence staining of NAIF1 in A431 cells. Cells were fixed with 4% PFA, permeabilized with 0.3% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-Human NAIF1 polyclonal antibody (dilution ratio 1:200) at 4 degrees celsius overnight. Then cells were stained with the Alexa Fluor (TM) 488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to Nucleus and Cytoplasm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-06197-100ul**

---

|                   |                                                     |
|-------------------|-----------------------------------------------------|
| HAF008            | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160            | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891        | Rabbit IgG Isotype Control                          |
| NBP2-51712-0.02mg | Recombinant Human NAIF1 His Protein                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-06197](http://www.novusbio.com/reviews/submit/NBP3-06197)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

