

# Product Datasheet

## CD200R1L Antibody (07) [HRP]

### NBP3-06393H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-06393H](http://www.novusbio.com/NBP3-06393H)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-06393H](http://www.novusbio.com/reviews/destination/NBP3-06393H)



**NBP3-06393H**

CD200R1L Antibody (07) [HRP]

| Product Information         |                                                                                                                                                                                                     |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                              |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                         |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                            |
| Clonality                   | Monoclonal                                                                                                                                                                                          |
| Clone                       | 07                                                                                                                                                                                                  |
| Preservative                | No Preservative                                                                                                                                                                                     |
| Isotype                     | IgG1                                                                                                                                                                                                |
| Conjugate                   | HRP                                                                                                                                                                                                 |
| Purity                      | Protein A purified                                                                                                                                                                                  |
| Buffer                      | PBS                                                                                                                                                                                                 |
| Product Description         |                                                                                                                                                                                                     |
| Host                        | Mouse                                                                                                                                                                                               |
| Gene ID                     | 344807                                                                                                                                                                                              |
| Gene Symbol                 | CD200R1L                                                                                                                                                                                            |
| Species                     | Cynomolgus Monkey                                                                                                                                                                                   |
| Reactivity Notes            | No cross-reactivity in ELISA with:<br/><br/>Human cell lysate (293 cell line)                                                                                                                       |
| Specificity/Sensitivity     | No cross-reactivity in ELISA with:<br><br>Human cell lysate (293 cell line)                                                                                                                         |
| Immunogen                   | This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Cynomolgus CD200R1L (Met1-Leu240). |
| Product Application Details |                                                                                                                                                                                                     |
| Applications                | ELISA                                                                                                                                                                                               |
| Recommended Dilutions       | ELISA                                                                                                                                                                                               |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                              |



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-06393H**

---

|                |                                        |
|----------------|----------------------------------------|
| NBP1-97074     | Mouse IgG1 Isotype Control (MG1) [HRP] |
| 7954-GM-010/CF | GM-CSF [Unconjugated]                  |
| 10198-CD-050   | CD200R1L [Unconjugated]                |
| NB100-39002    | FoxP3 Antibody - BSA Free              |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-06393H](http://www.novusbio.com/reviews/submit/NBP3-06393H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

