

# Product Datasheet

## Von Willebrand Factor Antibody (05) - Azide and BSA Free NBP3-06542-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-06542](http://www.novusbio.com/NBP3-06542)

Updated 2/26/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-06542](http://www.novusbio.com/reviews/destination/NBP3-06542)



**NBP3-06542-100ul**

Von Willebrand Factor Antibody (05) - Azide and BSA Free

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 100 ul                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 05                                                                                          |
| Preservative        | No Preservative                                                                             |
| Isotype             | IgG                                                                                         |
| Purity              | Protein A purified                                                                          |
| Buffer              | 0.2 um filtered solution in PBS                                                             |

| Product Description |                                                                                                                                                                                                                                      |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | This antibody can be stored at 2C-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Avoid repeated freeze-thaw cycles.          |
| Host                | Mouse                                                                                                                                                                                                                                |
| Gene ID             | 7450                                                                                                                                                                                                                                 |
| Gene Symbol         | VWF                                                                                                                                                                                                                                  |
| Species             | Human                                                                                                                                                                                                                                |
| Immunogen           | This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Von Willebrand Factor (Accession#: AAB59458.1; Met1-Lys2813). |

| Product Application Details |                                                                                                                                                                                                      |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | ELISA, Sandwich ELISA Capture                                                                                                                                                                        |
| Recommended Dilutions       | ELISA, Sandwich ELISA Capture                                                                                                                                                                        |
| Application Notes           | This antibody will detect Human Von Willebrand Factor in ELISA pair set (Catalog: # NBP2-79292). In a sandwich ELISA, it can be used as a capture antibody when paired with (Catalog: # NBP2-89724). |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-06542-100ul**

---

|                |                                                         |
|----------------|---------------------------------------------------------|
| HAF007         | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B        | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97019-5mg | Mouse IgG Isotype Control                               |
| NBP1-84762PEP  | Von Willebrand Factor Recombinant Protein Antigen       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-06542](http://www.novusbio.com/reviews/submit/NBP3-06542)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

