

# Product Datasheet

## Bovine Haptoglobin ELISA Kit (Colorimetric) NBP3-07118

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-07118](http://www.novusbio.com/NBP3-07118)

Updated 10/25/2022 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-07118](http://www.novusbio.com/reviews/destination/NBP3-07118)



**NBP3-07118****Bovine Haptoglobin ELISA Kit (Colorimetric)****Product Information**

<b>Unit Size</b>	1 Kit
<b>Concentration</b>	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
<b>Storage</b>	Storage of components varies. See protocol for specific instructions.

**Product Description**

<b>Gene ID</b>	3240
<b>Gene Symbol</b>	HP
<b>Species</b>	Bovine
<b>Kit Components</b>	Bovine Haptoglobin Microplate (1 EA), Sealing Tapes (3 EA), Bovine Haptoglobin Standard (60 ng lyophilized), Biotinylated Bovine Haptoglobin Antibody - 50X (120 uL), EIA Diluent Concentrate - 10X (30 mL), Wash Buffer Concentrate - 20X (30 mL), SP Conjugate - 100X (80 uL), Chromogen Substrate - 1X (8 mL), Stop Solution - 1X (12 mL)
<b>Standard Curve Range</b>	0.781 - 50 ng/mL (example only; lot dependent)
<b>Sensitivity</b>	0.65 ng/mL (example only; lot dependent)
<b>Inter Assay</b>	%CV < 5.6 (example only; lot dependent)
<b>Intra Assay</b>	%CV < 10.0 (example only; lot dependent)
<b>Assay Type</b>	Sandwich ELISA
<b>Spiking Recovery</b>	97% (example only; lot dependent)
<b>Suitable Sample Type</b>	Serum, Plasma, Cell Culture, Urine


**Product Application Details**

<b>Applications</b>	ELISA, Sandwich ELISA
<b>Recommended Dilutions</b>	ELISA, Sandwich ELISA

**Images**

ELISA: Bovine Haptoglobin ELISA Kit (Colorimetric) [NBP3-07118]

Bovine Haptoglobin Standard Curve





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-07118**

---

NBP2-51611-0.05mg	Recombinant Human Haptoglobin His Protein
210-TA-005	TNF-alpha [Unconjugated]
DHAPG0	Haptoglobin [HRP]
D6050	IL-6 [HRP]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-07118](http://www.novusbio.com/reviews/submit/NBP3-07118)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

