

# Product Datasheet

## CD79B Antibody (B29/123) [HRP] NBP3-08356H

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-08356H](http://www.novusbio.com/NBP3-08356H)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-08356H](http://www.novusbio.com/reviews/destination/NBP3-08356H)



**NBP3-08356H**

CD79B Antibody (B29/123) [HRP]

| Product Information         |                                                                                                                   |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                            |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                       |
| Storage                     | Store at 4C in the dark.                                                                                          |
| Clonality                   | Monoclonal                                                                                                        |
| Clone                       | B29/123                                                                                                           |
| Preservative                | No Preservative                                                                                                   |
| Isotype                     | IgG2b Kappa                                                                                                       |
| Conjugate                   | HRP                                                                                                               |
| Purity                      | Protein A or G purified                                                                                           |
| Buffer                      | PBS                                                                                                               |
| Product Description         |                                                                                                                   |
| Host                        | Mouse                                                                                                             |
| Gene ID                     | 974                                                                                                               |
| Gene Symbol                 | CD79B                                                                                                             |
| Species                     | Human                                                                                                             |
| Marker                      | B-Cell Marker                                                                                                     |
| Immunogen                   | Synthetic peptide representing the c-terminus from residue 215 of the murine CD79B polypeptide. (Uniprot: P40259) |
| Product Application Details |                                                                                                                   |
| Applications                | Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin         |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-08356H**

---

|                 |                                                            |
|-----------------|------------------------------------------------------------|
| NBP1-43317H     | Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) [HRP] |
| NBP2-59563-50ug | Recombinant Human CD79B His Protein                        |
| 7268-CT-100     | CTLA-4 [Unconjugated]                                      |
| 9687-CD-050     | CD79B [Unconjugated]                                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-08356H](http://www.novusbio.com/reviews/submit/NBP3-08356H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

