

Product Datasheet

ZAP70 Antibody (ZAP70/4040R) - Azide and BSA Free NBP3-08968

Unit Size: 100 ug

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-08968

Updated 7/16/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-08968



NBP3-08968

ZAP70 Antibody (ZAP70/4040R) - Azide and BSA Free

| Product Information | |
|-------------------------|---|
| Unit Size | 100 ug |
| Concentration | 1 mg/ml |
| Storage | Store at -20 to -80C. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | ZAP70/4040R |
| Preservative | No Preservative |
| Isotype | IgG |
| Purity | Protein A or G purified |
| Buffer | 1 mM PBS |
| Target Molecular Weight | 70 kDa |


| Product Description | |
|---------------------|---|
| Description | 1.0 mg/ml of antibody purified from Bioreactor Concentrate by Protein A/G. Prepared in 1mM PBS WITHOUT BSA & azide. Also available at 200 ug/ml WITH BSA & azide (NBP3-07223). Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80 C. |
| Host | Rabbit |
| Gene ID | 7535 |
| Gene Symbol | ZAP70 |
| Species | Human |
| Marker | Chronic Lymphocytic Leukemia Marker |
| Immunogen | Synthetic peptide corresponding to residues within aa200-350 of ZAP70 was used as an immunogen. (Uniprot: P43403) |

| Product Application Details | |
|-----------------------------|--|
| Applications | ELISA, Immunohistochemistry-Paraffin |
| Recommended Dilutions | ELISA 2-4 ug/ml, Immunohistochemistry-Paraffin 1-2 ug/ml |
| Application Notes | ELISA: Use antibody at 2-4ug/ml. For coating, order antibody without BSA |




Images

Immunohistochemistry-Paraffin: ZAP70 Antibody (ZAP70/4040R) - Azide and BSA Free [NBP3-08968] - Formalin-fixed, paraffin-embedded human tonsil stained with ZAP70 Recombinant Rabbit Monoclonal Antibody (ZAP70/4040R).



Immunohistochemistry-Paraffin: ZAP70 Antibody (ZAP70/4040R) - Azide and BSA Free [NBP3-08968] - Formalin-fixed, paraffin-embedded human tonsil stained with ZAP70 Recombinant Rabbit Monoclonal Antibody (ZAP70/4040R).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-08968

| | |
|---------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP1-87000PEP | ZAP70 Recombinant Protein Antigen |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-08968

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

