

# Product Datasheet

## TIAL1 Antibody (6E3) - BSA Free NBP3-11834

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-11834](http://www.novusbio.com/NBP3-11834)

Updated 5/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-11834](http://www.novusbio.com/reviews/destination/NBP3-11834)



**NBP3-11834**

TIAL1 Antibody (6E3) - BSA Free

**Product Information**

|                      |                             |
|----------------------|-----------------------------|
| <b>Unit Size</b>     | 0.1 mg                      |
| <b>Concentration</b> | 1 mg/ml                     |
| <b>Storage</b>       | Store at 4C. Do not freeze. |
| <b>Clonality</b>     | Monoclonal                  |
| <b>Clone</b>         | 6E3                         |
| <b>Preservative</b>  | 0.9% Sodium Azide           |
| <b>Isotype</b>       | IgG2a Kappa                 |
| <b>Purity</b>        | Protein A purified          |
| <b>Buffer</b>        | PBS                         |

**Product Description**


|                                |                                                                                                                                       |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>                    | Mouse                                                                                                                                 |
| <b>Gene ID</b>                 | 7073                                                                                                                                  |
| <b>Gene Symbol</b>             | TIAL1                                                                                                                                 |
| <b>Species</b>                 | Human, Mouse, Primate                                                                                                                 |
| <b>Specificity/Sensitivity</b> | The mouse monoclonal antibody 6E3 recognizes TIAL1, an ubiquitously expressed approximately 42 kDa intracellular RNA-binding protein. |
| <b>Immunogen</b>               | GST-tagged recombinant human TIAL1                                                                                                    |

**Product Application Details**

|                              |                                                                                                                                                   |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot, ELISA, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin                 |
| <b>Recommended Dilutions</b> | Western Blot, Flow Cytometry 5 ug/ml, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Flow (Intracellular) 5 ug/ml |

**Images**

Flow Cytometry: TIAL1 Antibody (6E3) [NBP3-11834] - Staining of TIAR in human cell line A-431 using purified mouse monoclonal antibody 6E3 (concentration in sample 5 ug/ml, GAM FITC, red-filled) from A-431 cells unstained by primary antibody (GAM FITC, black-dashed) in flow cytometry analysis (intracellular staining).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-11834**

|                  |                                                         |
|------------------|---------------------------------------------------------|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B          | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-96981-0.5mg | Mouse IgG2a Kappa Isotype Control (M2AK)                |
| NBP2-51914-0.1mg | Recombinant Human TIAL1 His Protein                     |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-11834](http://www.novusbio.com/reviews/submit/NBP3-11834)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

