

# Product Datasheet

## PGAP3 Antibody NBP3-12195

Unit Size: 100 ug

Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-12195](http://www.novusbio.com/NBP3-12195)

Updated 1/2/2022 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-12195](http://www.novusbio.com/reviews/destination/NBP3-12195)



**NBP3-12195****PGAP3 Antibody**

| <b>Product Information</b> |                                                                                                                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 100 ug                                                                                                                                                                                    |
| <b>Concentration</b>       | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                               |
| <b>Storage</b>             | Store at 4C short term. Store at -20C long term. Avoid freeze-thaw cycles.                                                                                                                |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                                                                                                        |
| <b>Purity</b>              | Affinity purified                                                                                                                                                                         |
| <b>Buffer</b>              | Contains Tris, HCl/glycine buffer (pH 7.4-7.8), 30% glycerol, 0.5% BSA, along with cryo-protective agents, and HEPES (exact storage and stabilization buffer information is proprietary). |

| <b>Product Description</b> |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Host</b>                | Rabbit                                                                         |
| <b>Gene ID</b>             | 93210                                                                          |
| <b>Gene Symbol</b>         | PGAP3                                                                          |
| <b>Species</b>             | Human, Primate                                                                 |
| <b>Immunogen</b>           | Synthetic peptide taken within amino acid region 10-60 on human PGAP3 protein. |

| <b>Product Application Details</b> |                                                                          |
|------------------------------------|--------------------------------------------------------------------------|
| <b>Applications</b>                | Western Blot, Dot Blot, ELISA, Immunohistochemistry, Immunoprecipitation |
| <b>Recommended Dilutions</b>       | Western Blot, ELISA, Immunohistochemistry, Immunoprecipitation, Dot Blot |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-12195**

---

|                   |                                                                        |
|-------------------|------------------------------------------------------------------------|
| HAF008            | Goat anti-Rabbit IgG Secondary Antibody [HRP (Horseradish Peroxidase)] |
| NB7160            | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]                    |
| H00093210-P01-2ug | Recombinant Human PGAP3 GST (N-Term) Protein                           |
| 664-LI-025/CF     | LIGHT/TNFSF14 [Unconjugated]                                           |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-12195](http://www.novusbio.com/reviews/submit/NBP3-12195)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

