

# Product Datasheet

## **RBM19 Antibody** **NBP3-12591-100ul**

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-12591](http://www.novusbio.com/NBP3-12591)

Updated 9/5/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-12591](http://www.novusbio.com/reviews/destination/NBP3-12591)



**NBP3-12591-100ul**


RBM19 Antibody

| Product Information         |                                                                                                                                                                                                                                                                                    |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                                                                                                                                                             |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                        |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                                                                                                                                                                             |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                         |
| Preservative                | 0.03% Proclin 300                                                                                                                                                                                                                                                                  |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                |
| Purity                      | Antigen and protein A Affinity-purified                                                                                                                                                                                                                                            |
| Buffer                      | PBS, pH7.0                                                                                                                                                                                                                                                                         |
| Product Description         |                                                                                                                                                                                                                                                                                    |
| Description                 | <p>This antibody can be stored at 2-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Avoid repeated freeze-thaw cycles.</p> <p>Protein A &amp; Antigen Affinity purified</p> |
| Host                        | Rabbit                                                                                                                                                                                                                                                                             |
| Gene Symbol                 | RBM19                                                                                                                                                                                                                                                                              |
| Species                     | Human                                                                                                                                                                                                                                                                              |
| Immunogen                   | Produced in rabbits immunized with E. coli-derived Human RBM19 fragment                                                                                                                                                                                                            |
| Product Application Details |                                                                                                                                                                                                                                                                                    |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunoprecipitation                                                                                                                                                                                                                   |
| Recommended Dilutions       | Western Blot 1:500-1:2000, Immunoprecipitation 1-5uL/mg of lysate, Immunohistochemistry-Paraffin 1:50-1:200                                                                                                                                                                        |




## Images


Western Blot: RBM19 Antibody [NBP3-12591] - Lane A: U-251 MG Whole Cell Lysate Lane B: 293T Whole Cell Lysate Lysates/proteins at 30 ug per lane. Secondary Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution. Developed using the ECL technique. Performed under reducing conditions. Predicted band size:107 kDa Observed band size:130 kDa



Immunohistochemistry-Paraffin: RBM19 Antibody [NBP3-12591] - Staining of human RBM19 in human brain with rabbit polyclonal antibody at 1:100 dilution.



Immunoprecipitation: RBM19 Antibody [NBP3-12591] - Lane A:0.5 mg U-251 MG Whole Cell Lysate 4 uL anti-RBM19 rabbit polyclonal antibody and 60 ug of Immunomagnetic beads Protein A/G. Primary antibody: Anti-RBM19 rabbit polyclonal antibody, at 1:100 dilution Secondary antibody: Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution Developed using the ECL technique. Performed under reducing conditions. Predicted band size: 107 kDa Observed band size :130 kDa



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-12591-100ul**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP1-84135PEP | RBM19 Recombinant Protein Antigen                   |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-12591](http://www.novusbio.com/reviews/submit/NBP3-12591)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

