

# Product Datasheet

## Influenza A H3N2 Hemagglutinin Antibody - (A/Brisbane/10/2007) NBP3-12797-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-12797](http://www.novusbio.com/NBP3-12797)

Updated 6/18/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-12797](http://www.novusbio.com/reviews/destination/NBP3-12797)



**NBP3-12797-100ul**

Influenza A H3N2 Hemagglutinin Antibody - (A/Brisbane/10/2007)

**Product Information**

|                      |                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 100 ul                                                                                      |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| <b>Clonality</b>     | Polyclonal                                                                                  |
| <b>Preservative</b>  | No Preservative                                                                             |
| <b>Isotype</b>       | IgG                                                                                         |
| <b>Purity</b>        | Protein A purified                                                                          |
| <b>Buffer</b>        | 0.2 um filtered solution in PBS                                                             |

**Product Description**


|                                |                                                                                                                                                                                                                                               |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>             | This antibody can be stored at 2-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Preservative-Free. Avoid repeated freeze-thaw cycles. |
| <b>Host</b>                    | Rabbit                                                                                                                                                                                                                                        |
| <b>Species</b>                 | Influenza A Virus H3N2                                                                                                                                                                                                                        |
| <b>Specificity/Sensitivity</b> | Influenza A H3N2 Hemagglutinin (A/Brisbane/10/2007)                                                                                                                                                                                           |
| <b>Immunogen</b>               | Purified, recombinant Influenza A H3N2 Hemagglutinin (ABW23353.1; Met 1-Ile 531)                                                                                                                                                              |

**Product Application Details**

|                              |                                                  |
|------------------------------|--------------------------------------------------|
| <b>Applications</b>          | Western Blot, ELISA                              |
| <b>Recommended Dilutions</b> | Western Blot 1:1000-1:5000, ELISA 1:5000-1:10000 |

**Images**

Western Blot: Influenza A H3N2 Hemagglutinin Antibody - (A/Brisbane/10/2007) [NBP3-12797] - Influenza A H3N2 Hemagglutinin (A/Brisbane/10/2007) Recombinant Protein Lane A: 30ng Lane B: 5ng Secondary Goat Anti- Rabbit IgG H&L (Dylight 800) at 1/10000 dilution. Developed using the ECL technique. Performed under reducing conditions.



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

**Products Related to NBP3-12797-100ul**

---

|                  |                                                                                       |
|------------------|---------------------------------------------------------------------------------------|
| HAF008           | Goat anti-Rabbit IgG Secondary Antibody [HRP]                                         |
| NB7156           | Goat anti-Rabbit IgG (H+L) Secondary Antibody                                         |
| NBP2-24891       | Rabbit IgG Isotype Control                                                            |
| NBP1-99037-100ug | Recombinant Influenza A H3N2 Hemagglutinin<br>(A/canine/Guangdong/1/2006) His Protein |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-12797](http://www.novusbio.com/reviews/submit/NBP3-12797)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

