

Product Datasheet

Influenza A H3N2 Hemagglutinin Antibody - (A/Aichi/2/1968) NBP3-12805-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-12805

Updated 6/18/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-12805



NBP3-12805-100ul

Influenza A H3N2 Hemagglutinin Antibody - (A/Aichi/2/1968)

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	0.2 um filtered solution in PBS

Product Description


Description	This antibody can be stored at 2-8C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20C to -80C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Species	Influenza A Virus H3N2
Reactivity Notes	Influenza A H3N2 (A/Aichi/2/1968)
Specificity/Sensitivity	Influenza A H3N2 Hemagglutinin (A/Aichi/2/1968)
Immunogen	Purified, recombinant Influenza A H3N2 Hemagglutinin (AAA43178.1; Met1-Arg345)

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 1:1000-1:5000

Images

Western Blot: Influenza A H3N2 Hemagglutinin Antibody - (A/Aichi/2/1968) [NBP3-12805] - Sample: Influenza A H3N2 Hemagglutinin (A/Aichi/2/1968) Recombinant Protein Lane A: 50ng Lane B: 10ng Secondary Goat Anti- Rabbit IgG H&L (Dylight 800) at 1/10000 dilution. Developed using the Odyssey technique. Performed under reducing conditions.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-12805-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7156	Goat anti-Rabbit IgG (H+L) Secondary Antibody
NBP2-24891	Rabbit IgG Isotype Control
NBP1-99037-100ug	Recombinant Influenza A H3N2 Hemagglutinin (A/canine/Guangdong/1/2006) His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-12805

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

