

Product Datasheet

LAG-3 Antibody (1C6) - VHH - Azide and BSA Free NBP3-12839

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-12839

Updated 7/31/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-12839



NBP3-12839

LAG-3 Antibody (1C6) - VHH - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1C6
Preservative	No Preservative
Purity	Affinity purified
Buffer	PBS


Product Description	
Description	Recombinant Monoclonal Antibody - Llamabody(TM). Llama VHH domain containing a Myc-tag for detection and His-tag (MW ~ 15kDa)
Host	Llama
Gene ID	3902
Gene Symbol	LAG3
Species	Human
Immunogen	LAG-3 antibody was raised against a recombinant protein corresponding to amino acids 23-434 of human LAG-3.

Product Application Details	
Applications	Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Flow Cytometry 1 ug/mL, ELISA 0.1 ug/mL, Immunohistochemistry 1 ug/mL, Immunocytochemistry/ Immunofluorescence 5 ug/mL, Immunohistochemistry-Paraffin




Images


Immunocytochemistry/Immunofluorescence: LAG-3 Antibody (1C6) - VHH [NBP3-12839] - Immunofluorescence of LAG-3 in transfected HEK293 cells with LAG-3 single domain antibody at 5 ug/mL.




Immunohistochemistry: LAG-3 Antibody (1C6) - VHH [NBP3-12839] - Immunohistochemistry of LAG-3 in human lymph node tissue with LAG-3 single domain antibody at 1 ug/mL.



Flow Cytometry: LAG-3 Antibody (1C6) - VHH [NBP3-12839] - Flow cytometry analysis of LAG-3 overexpressing HEK293 cells using LAG-3 single domain antibody at 1 ug/mL. Blue: Untransfected HEK293 cells. Yellow: LAG-3 overexpressing HEK293 cells.



ELISA: LAG-3 Antibody (1C6) - VHH [NBP3-12839] - Titration ELISA analysis of LAG-3 sdAbs to detect recombinant LAG-3 (extracellular domain) coated at 1 ug/mL. sdAbs are detected with a mouse mAb against a C-terminal myc-tag followed by a goat anti-mouse IgG-HRP conjugate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-12839

AF011-500	Goat anti-Llama IgG Secondary Antibody [Unconjugated]
NB7240	Goat anti- Alpaca, Llama IgG (H+L) Secondary Antibody
NBP1-85781PEP	LAG-3 Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-12839

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

