

Product Datasheet

UQCRB Antibody - BSA Free

NBP3-12970

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-12970

Updated 3/3/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-12970



NBP3-12970

UQCRB Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol

Product Description

Description	Centrifuge briefly prior to opening.
Host	Rabbit
Gene ID	7381
Gene Symbol	UQCRB
Species	Human, Mouse, Rat
Reactivity Notes	Mouse (88%), Pig (88%), Chicken (81%), Bovine (86%)
Immunogen	Recombinant protein encompassing a sequence within the center region of human UQCRB. The exact sequence is proprietary.


Product Application Details

Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin




Images


Western Blot: UQCRB Antibody [NBP3-12970] - UQCRB antibody detects UQCRB protein by western blot analysis. Human tissue extracts (30 ug) and mouse tissue extracts (50 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with UQCRB antibody (NBP3-12970) diluted at 1:5000.



Immunohistochemistry-Paraffin: UQCRB Antibody [NBP3-12970] - UQCRB antibody detects UQCRB protein at mitochondria in rat intestine by immunohistochemical analysis. Sample: Paraffin-embedded rat intestine. UQCRB antibody (NBP3-12970) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: UQCRB Antibody [NBP3-12970] - UQCRB antibody detects UQCRB protein at mitochondria in mouse liver by immunohistochemical analysis. Sample: Paraffin-embedded mouse liver. UQCRB antibody (NBP3-12970) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-12970

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00007381-P01-10ug	Recombinant Human UQCRB GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-12970

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

