

# Product Datasheet

## TECPR1 Antibody - BSA Free

### NBP3-13113

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-13113](http://www.novusbio.com/NBP3-13113)

Updated 3/3/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-13113](http://www.novusbio.com/reviews/destination/NBP3-13113)



**NBP3-13113**

TECPR1 Antibody - BSA Free

**Product Information**

|                                |                                                                                                                  |
|--------------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>               | 100 ul                                                                                                           |
| <b>Concentration</b>           | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>                 | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| <b>Clonality</b>               | Polyclonal                                                                                                       |
| <b>Preservative</b>            | 0.025% Proclin 300                                                                                               |
| <b>Isotype</b>                 | IgG                                                                                                              |
| <b>Purity</b>                  | Protein A purified                                                                                               |
| <b>Buffer</b>                  | PBS (pH 7), 20% Glycerol                                                                                         |
| <b>Target Molecular Weight</b> | 130 kDa                                                                                                          |

**Product Description**


|                    |                                                                                                                              |
|--------------------|------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Centrifuge briefly prior to opening.                                                                                         |
| <b>Host</b>        | Rabbit                                                                                                                       |
| <b>Gene ID</b>     | 25851                                                                                                                        |
| <b>Gene Symbol</b> | TECPR1                                                                                                                       |
| <b>Species</b>     | Mouse                                                                                                                        |
| <b>Immunogen</b>   | Recombinant protein encompassing a sequence within the N-terminus region of human TECPR1. The exact sequence is proprietary. |

**Product Application Details**

|                              |                           |
|------------------------------|---------------------------|
| <b>Applications</b>          | Western Blot              |
| <b>Recommended Dilutions</b> | Western Blot 1:500-1:3000 |

**Images**

Western Blot: TECPR1 Antibody [NBP3-13113] - Various whole cell and tissue extracts (60 ug) were separated by 5% SDS-PAGE, and the membrane was blotted with TECPR1 antibody (NBP3-13113) diluted at 1:1000.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-13113**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP1-84206PEP | TECPR1 Recombinant Protein Antigen                  |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-13113](http://www.novusbio.com/reviews/submit/NBP3-13113)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

