

Product Datasheet

Methoxypolyethylene Glycol Antibody (HL1020) - Azide and BSA Free NBP3-13681

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13681

Updated 7/31/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13681



NBP3-13681


Methoxypolyethylene Glycol Antibody (HL1020) - Azide and BSA Free

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL1020
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS
Product Description	
Description	Centrifuge briefly prior to opening.
Host	Rabbit
Species	Non-species specific
Immunogen	The immunogen used to generate this antibody corresponds to Methoxypolyethylene Glycol (methoxy group).
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:5000-1:128000, ELISA Assay dependent




Images

Western Blot: Methoxypolyethylene Glycol Antibody (HL1020) [NBP3-13681] - Polyethylene Glycolylated OVA (1 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with Polyethylene Glycol (methoxy group) antibody [HL1020] (NBP3-13681) diluted at 1:128000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



ELISA: Methoxypolyethylene Glycol Antibody (HL1020) [NBP3-13681] - Indirect ELISA analysis performed by coating plate with recombinant OVA, Polyethylene Glycolylated OVA (Polyethylene Glycol24), Polyethylene Glycolylated OVA (Polyethylene Glycol3), and Polyethylene Glycolylated OVA acid (Polyethylene Glycol3 acid) (100-1.56 ng/mL). Coated protein was probed with Polyethylene Glycol (methoxy group) antibody [HL1020] (NBP3-13681) (1 ug/mL). Rabbit IgG antibody (HRP) (NBP2-19301) (1:10000) was used to detect bound primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-13681

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13681

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

