

Product Datasheet

RBPMS Antibody (HL1106) - Azide and BSA Free NBP3-13723

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13723

Updated 7/31/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13723



NBP3-13723

RBPMS Antibody (HL1106) - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HL1106
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS

Product Description


Description	Centrifuge briefly prior to opening.
Host	Rabbit
Gene ID	11030
Gene Symbol	RBPMS
Species	Human, Mouse, Rat
Immunogen	Recombinant protein encompassing a sequence within the center region of human RBPMS. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot
Recommended Dilutions	Western Blot

Images


Western Blot: RBPMS Antibody (HL1106) [NBP3-13723] - Various tissue extracts (50 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with RBPMS antibody [HL1106] (NBP3-13723) diluted at 1:1000. The HRP-conjugated anti-FAHSE was used to detect the primary antibody.




Western Blot: RBPMS Antibody (HL1106) [NBP3-13723] - Various whole cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with RBPMS antibody [HL1106] (NBP3-13723) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Immunohistochemistry-Paraffin: RBPMS Antibody (HL1106) [NBP3-13723] - RBPMS antibody [HL1106] detects RBPMS protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded mouse intestine. RBPMS stained by RBPMS antibody [HL1106] (NBP3-13723) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




Immunohistochemistry-Paraffin: RBPMS Antibody (HL1106) [NBP3-13723] - RBPMS antibody [HL1106] detects RBPMS protein by immunohistochemical analysis. Sample: Paraffin-embedded mouse tissues. RBPMS stained by RBPMS antibody [HL1106] (NBP3-13723) diluted at 1:100. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




Immunohistochemistry-Paraffin: RBPMS Antibody (HL1106) [NBP3-13723] - RBPMS antibody [HL1106] detects RBPMS protein at cytoplasm and nucleus by immunohistochemical analysis. Sample: Paraffin-embedded mouse eye. Green: RBPMS stained by RBPMS antibody [HL1106] (NBP3-13723) diluted at 1:250. Red: beta Tubulin 3/ Tuj1, a neural marker, stained by beta Tubulin 3/ Tuj1 antibody [GT11710] diluted at 1:500. Blue: Fluoroshield with DAPI. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: RBPMS Antibody (HL1106) [NBP3-13723] - RBPMS antibody [HL1106] detects RBPMS protein at cytoplasm and nucleus by immunohistochemical analysis. Sample: Paraffin-embedded rat eye. Green: RBPMS stained by RBPMS antibody [HL1106] (NBP3-13723) diluted at 1:250. Red: beta Tubulin 3/ Tuj1, a neural marker, stained by beta Tubulin 3/ Tuj1 antibody [GT11710] diluted at 1:500. Blue: Fluoroshield with DAPI. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Western Blot: RBPMS Antibody (HL1106) [NBP3-13723] - Non-transfected (–) and transfected (+) HCT-116 whole cell extract (60 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with RBPMS antibody [HL1106] diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-13723

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00011030-P01-10ug	Recombinant Human RBPMS GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13723

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

