

Product Datasheet

IGF-I R/IGF1R Antibody (IGF1R/4667) NBP3-13948-100ug

Unit Size: 100 ug

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-13948

Updated 7/16/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-13948



NBP3-13948-100ug

IGF-I R/IGF1R Antibody (IGF1R/4667)

Product Information	
Unit Size	100 ug
Concentration	0.2 mg/ml
Storage	Store at 4C.
Clonality	Monoclonal
Clone	IGF1R/4667
Preservative	0.05% Sodium Azide
Isotype	IgG2c Kappa
Purity	Protein A or G purified
Buffer	10 mM PBS with 0.05% BSA


Product Description	
Description	<p>200ug/ml of antibody purified from Bioreactor Concentrate by Protein A or G. Prepared in 10 mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0 mg/ml. (NBP3-14094)</p> <p>Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80 C. Non-hazardous.</p> <p>Theoretical molecular weight: 130kDa (IGF-1R); 97kDa (IGF-1R); 200kDa (pro-IGF-1R)</p>
Host	Mouse
Gene ID	3480
Gene Symbol	IGF1R
Species	Human
Immunogen	Recombinant fragment (around aa550-750) of human IGF-I R/IGF1R protein (exact sequence is proprietary) (Uniprot: P08069)

Product Application Details	
Applications	Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, In vitro assay, Protein Array
Recommended Dilutions	Flow Cytometry 1-2 ug/million cells, Immunocytochemistry/ Immunofluorescence 1-4 ug/ml, Immunohistochemistry-Paraffin, In vitro assay, Protein Array
Application Notes	In-vitro Neutralization of Biological Activity of IGF-1: Order Ab without Azide




Images

Immunohistochemistry-Paraffin: IGF-I R/IGF1R Antibody (IGF1R/4667) [NBP3-13948] - Formalin-fixed, paraffin-embedded human prostate carcinoma stained with IGF-I R/IGF1R antibody (IGF1R/4667).




Immunohistochemistry-Paraffin: IGF-I R/IGF1R Antibody (IGF1R/4667) [NBP3-13948] - Formalin-fixed, paraffin-embedded human colon carcinoma stained with IGF-I R/IGF1R antibody (IGF1R/4667).



Immunohistochemistry-Paraffin: IGF-I R/IGF1R Antibody (IGF1R/4667) [NBP3-13948] - Formalin-fixed, paraffin-embedded human breast carcinoma stained with IGF-I R/IGF1R antibody (IGF1R/4667).



Protein Array: IGF-I R/IGF1R Antibody (IGF1R/4667) [NBP3-13948] - Analysis of Protein Array containing more than 19,000 full-length human proteins using IGF-I R/IGF1R Antibody (IGF1R/4667).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-13948-100ug

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-77680PEP	IGF-I R/IGF1R Antibody Blocking Peptide
236-EG-200	EGF [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-13948

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

