

# Product Datasheet

## ZC3H7A Antibody (PCRP-ZC3H7A-1D6) [DyLight 594] NBP3-14047DL594

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14047DL594](http://www.novusbio.com/NBP3-14047DL594)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14047DL594](http://www.novusbio.com/reviews/destination/NBP3-14047DL594)



**NBP3-14047DL594**

ZC3H7A Antibody (PCRP-ZC3H7A-1D6) [DyLight 594]

| Product Information         |                                                                                                   |
|-----------------------------|---------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                            |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.       |
| Storage                     | Store at 4C in the dark.                                                                          |
| Clonality                   | Monoclonal                                                                                        |
| Clone                       | PCRP-ZC3H7A-1D6                                                                                   |
| Preservative                | 0.05% Sodium Azide                                                                                |
| Isotype                     | IgG1                                                                                              |
| Conjugate                   | DyLight 594                                                                                       |
| Purity                      | Protein A or G purified                                                                           |
| Buffer                      | 50mM Sodium Borate                                                                                |
| Product Description         |                                                                                                   |
| Host                        | Mouse                                                                                             |
| Gene ID                     | 29066                                                                                             |
| Gene Symbol                 | ZC3H7A                                                                                            |
| Species                     | Human                                                                                             |
| Immunogen                   | Recombinant full-length human ZC3H7A protein (Uniprot: Q8IWR0)                                    |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                 |
| Product Application Details |                                                                                                   |
| Applications                | ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Protein Array |
| Recommended Dilutions       | Flow Cytometry, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Protein Array |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-14047DL594**

---

|                    |                                                  |
|--------------------|--------------------------------------------------|
| NBP1-97005DL594    | Mouse IgG1 Isotype Control (MG1) [DyLight 594]   |
| H00029066-P01-10ug | Recombinant Human ZC3H7A GST (N-Term) Protein    |
| H00029066-T01      | ZC3H7A 293T Cell Transient Overexpression Lysate |
| NB200-103          | p53 Antibody (PAb 240) - BSA Free                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14047DL594](http://www.novusbio.com/reviews/submit/NBP3-14047DL594)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

