

Product Datasheet

DNA Polymerase epsilon subunit 3 Antibody (PCRP-POLE3-3D3) [mFluor Violet 450 SE] NBP3-14143MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-14143MFV450

Updated 10/26/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-14143MFV450



NBP3-14143MFV450

DNA Polymerase epsilon subunit 3 Antibody (PCRP-POLE3-3D3) [mFluor Violet 450 SE]


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	PCRP-POLE3-3D3
Preservative	0.05% Sodium Azide
Isotype	IgG2a
Conjugate	mFluor Violet 450 SE
Purity	Protein A or G purified
Buffer	50mM Sodium Borate

Product Description	
Host	Mouse
Gene ID	54107
Gene Symbol	POLE3
Species	Human
Reactivity Notes	Predicted to react in Mouse, Rat and Xenopus.
Immunogen	Recombinant full-length human DNA Polymerase epsilon subunit 3 protein (Uniprot: Q9NRF9)
Notes	mFluor(TM) is a trademark of AAT Bioquest, Inc. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.

Product Application Details	
Applications	ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Protein Array
Recommended Dilutions	Flow Cytometry, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Protein Array
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

DNA Polymerase epsilon subunit 3 Antibody (PCRP-POLE3-3D3) [mFluor Violet 450 SE] [NBP3-14143MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-14143MFV450

NBP2-23341	Recombinant Human DNA Polymerase epsilon subunit 3 His Protein
DRT100	TNF RI/TNFRSF1A [HRP]
H00054107-T01	DNA Polymerase epsilon subunit 3 293T Cell Transient Overexpression Lysate
AF835	Caspase-3 Antibody [Unconjugated] - Active

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-14143MFV450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



