

# Product Datasheet

## ZHX3 Antibody (PCRP-ZHX3-1D11) [PE] NBP3-14225PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14225PE](http://www.novusbio.com/NBP3-14225PE)

Updated 10/26/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14225PE](http://www.novusbio.com/reviews/destination/NBP3-14225PE)



**NBP3-14225PE**

ZHX3 Antibody (PCRP-ZHX3-1D11) [PE]

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.                                                                    |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | PCRP-ZHX3-1D11                                                                              |
| Preservative                | 0.05% Sodium Azide                                                                          |
| Isotype                     | IgG2b                                                                                       |
| Conjugate                   | PE                                                                                          |
| Purity                      | Protein A or G purified                                                                     |
| Buffer                      | PBS                                                                                         |
| Product Description         |                                                                                             |
| Host                        | Mouse                                                                                       |
| Gene ID                     | 23051                                                                                       |
| Gene Symbol                 | ZHX3                                                                                        |
| Species                     | Human                                                                                       |
| Reactivity Notes            | 0                                                                                           |
| Immunogen                   | Recombinant full-length human ZHX3 protein (Uniprot: Q9H4I2)                                |
| Product Application Details |                                                                                             |
| Applications                | Western Blot, ELISA, Flow Cytometry, Immunoprecipitation, Protein Array                     |
| Recommended Dilutions       | Western Blot, Flow Cytometry, ELISA, Immunoprecipitation, Protein Array                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                      |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-14225PE**

---

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| NBP2-27231PE       | Mouse IgG2b Isotype Control (MPC-11) [PE]              |
| H00023051-Q01-10ug | Recombinant Human ZHX3 GST (N-Term) Protein            |
| H00023051-T01      | ZHX3 293T Cell Transient Overexpression Lysate         |
| MAB1368            | alpha-Fetoprotein/AFP Antibody (189502) [Unconjugated] |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14225PE](http://www.novusbio.com/reviews/submit/NBP3-14225PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

