

# Product Datasheet

## Human Fas Ligand/TNFSF6 ELISA Kit (Colorimetric) NBP3-14629

Unit Size: 1 Kit

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14629](http://www.novusbio.com/NBP3-14629)

Updated 1/1/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14629](http://www.novusbio.com/reviews/destination/NBP3-14629)



**NBP3-14629****Human Fas Ligand/TNFSF6 ELISA Kit (Colorimetric)**


<b>Product Information</b>	
<b>Unit Size</b>	1 Kit
<b>Concentration</b>	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
<b>Storage</b>	Store at 4C. Do not freeze.

<b>Product Description</b>	
<b>Description</b>	Incubation: 3 hours 45 minutes
<b>Gene ID</b>	356
<b>Gene Symbol</b>	FASLG
<b>Species</b>	Human
<b>Specificity/Sensitivity</b>	Recognizes both natural and recombinant soluble human Fas Ligand/TNFSF6
<b>Kit Components</b>	Pre-coated strip plates, Biotinylated secondary antibody, Standards, Controls (where applicable), Buffers, Streptavidin-HRP, TMB, Stop reagents, Detailed procedure.
<b>Standard Curve Range</b>	62.5-2000 pg/ml
<b>Sensitivity</b>	12 pg/ml
<b>Intra Assay</b>	CV% < 4.5
<b>Assay Type</b>	Sandwich ELISA
<b>Suitable Sample Type</b>	Serum, Plasma, Cell culture supernatant
<b>Sample Volume</b>	100 ul

<b>Product Application Details</b>	
<b>Applications</b>	ELISA
<b>Recommended Dilutions</b>	ELISA

**Images**

ELISA: Human Fas Ligand/TNFSF6 ELISA Kit (Colorimetric) [NBP3-14629] - Example of Human Fas Ligand/TNFSF6 Standard Curve.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-14629**

---

NBP2-49160PEP	Fas Ligand/TNFSF6 Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
DFL00B	Fas Ligand/TNFSF6 [HRP]
285-IF-100	IFN-gamma [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14629](http://www.novusbio.com/reviews/submit/NBP3-14629)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

