

Product Datasheet

FoxP3 Antibody (BLR034F)

NBP3-14670

Unit Size: 100 ul

Store at 2 - 8 C / 1 year from date of receipt

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-14670

Updated 8/11/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-14670



NBP3-14670

FoxP3 Antibody (BLR034F)

Product Information

| | |
|----------------------|---|
| Unit Size | 100 ul |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 2 - 8 C / 1 year from date of receipt |
| Clonality | Monoclonal |
| Clone | BLR034F |
| Preservative | 0.09% Sodium Azide |
| Purity | >95% |
| Buffer | Borate Buffered Saline (BBS) pH 8.2 with , BSA Free |

Product Description


| | |
|--------------------|---------------------|
| Host | Rabbit |
| Gene ID | 50943 |
| Gene Symbol | FOXP3 |
| Species | Human |
| Immunogen | between 275 and 325 |

Product Application Details


| | |
|------------------------------|---|
| Applications | Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Multiplex Immunofluorescence |
| Recommended Dilutions | Western Blot 1:1000, Immunohistochemistry 1:100 - 1:500, Immunocytochemistry/ Immunofluorescence 1:100 - 1:500, Immunoprecipitation 20 ul/1mg lysate, Immunohistochemistry-Paraffin 1:100 - 1:500, Multiplex Immunofluorescence 1:250 |
| Application Notes | Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE cell sections. |

Images


Western Blot: FoxP3 Antibody (BLR034F) [NBP3-14670] - Whole cell lysate (50 ug) from Jurkat, HEK293T, MJ, MOLT-4, and RPMI-8226 cells prepared using NETN lysis buffer. Rabbit anti-FOXP3 recombinant monoclonal antibody (NBP3-14670 lot 2) used at 1:1000. Secondary: HRP-conjugated goat anti-rabbit IgG (NB7156). Chemiluminescence with an exposure time of 30 seconds. Lower Panel: Rabbit anti-COPB2 antibody.




Immunohistochemistry: FoxP3 Antibody (BLR034F) [NBP3-14670] - FFPE section of human breast carcinoma. Antibody: Rabbit anti-FOXP3 recombinant monoclonal antibody used at 1:250. Secondary: HRP-conjugated goat antirabbit IgG. Substrate: Opal™. Counterstain: DAPI (blue).




Immunohistochemistry-Paraffin: FoxP3 Antibody (BLR034F) [NBP3-14670] - Antibody: Rabbit anti-FOXP3 recombinant monoclonal antibody. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.




Immunohistochemistry-Paraffin: FoxP3 Antibody (BLR034F) [NBP3-14670] - Antibody: Rabbit anti-FOXP3 recombinant monoclonal antibody. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.



Immunohistochemistry: FoxP3 Antibody (BLR034F) [NBP3-14670] - Antibody: Rabbit anti-FOXP3 recombinant monoclonal antibody. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.



Immunoprecipitation: FoxP3 Antibody (BLR034F) [NBP3-14670] -
Samples: Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded)
from MJ cells prepared using NETN lysis buffer. Antibodies: Rabbit anti-
FOXP3 recombinant monoclonal antibody used for IP at 20 μ l/mg lysate.
FOXP3 was also immunoprecipitated by a previous lot of this antibody.
For blotting immunoprecipitated FOXP3, was used at 1:1000.
Chemiluminescence with an exposure time of 30 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-14670

| | |
|--------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NB600-246PEP | FoxP3 Antibody Blocking Peptide |
| 210-TA-005 | TNF-alpha [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-14670

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

