

# Product Datasheet

## Twist-1 Antibody NBP3-14732

Unit Size: 100 ul

Store at 2 - 8 C / 1 year from date of receipt

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-14732](http://www.novusbio.com/NBP3-14732)

Updated 4/20/2022 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-14732](http://www.novusbio.com/reviews/destination/NBP3-14732)



**NBP3-14732**

## Twist-1 Antibody

**Product Information**

<b>Unit Size</b>	100 ul
<b>Concentration</b>	Please see the vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at 2 - 8 C / 1 year from date of receipt
<b>Clonality</b>	Monoclonal
<b>Preservative</b>	0.09% Sodium Azide
<b>Purity</b>	>95%
<b>Buffer</b>	Phosphate Buffered Saline (PBS) containing 0.2% BSA and

**Product Description**


<b>Host</b>	Rabbit
<b>Gene ID</b>	7291
<b>Gene Symbol</b>	TWIST1
<b>Species</b>	Human
<b>Immunogen</b>	between 1 and 50

**Product Application Details**

<b>Applications</b>	Western Blot, Immunocytochemistry/Immunofluorescence, Immunoprecipitation
<b>Recommended Dilutions</b>	Western Blot, Immunocytochemistry/Immunofluorescence, Immunoprecipitation

**Images**

Immunoprecipitation: Twist-1 Antibody [NBP3-14732] - Whole cell lysate (1.0 mg per IP reaction; 20% of IP loaded) from SK-N-SH cells prepared using NETN lysis buffer. Antibodies: Rabbit antiTWIST1 recombinant monoclonal antibody used for IP at 20 ul/mg lysate. TWIST1 was also immunoprecipitated by a second antibody against a different epitope of TWIST1. For blotting immunoprecipitated TWIST1, was used at 1:1000. Chemiluminescence with an exposure time of 10 seconds.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-14732**

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB7156	Goat anti-Rabbit IgG (H+L) Secondary Antibody
H00007291-Q01-10ug	Recombinant Human Twist-1 GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-14732](http://www.novusbio.com/reviews/submit/NBP3-14732)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

