

Product Datasheet

PI 3-Kinase p85 alpha Antibody (S08-8H6) NBP3-14962-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-14962

Updated 4/15/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-14962



NBP3-14962-100ul

PI 3-Kinase p85 alpha Antibody (S08-8H6)

| Product Information | |
|-------------------------|---|
| Unit Size | 100 ul |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | S08-8H6 |
| Preservative | 0.01% Sodium Azide |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA |
| Target Molecular Weight | 84 kDa |

| Product Description | |
|---------------------|---|
| Host | Rabbit |
| Gene Symbol | PIK3R1 |
| Species | Human, Mouse, Rat |
| Immunogen | A synthetic peptide of human PI 3-Kinase p85 alpha (Uniprot#: P27986) |

| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions | Western Blot 1:1000, ELISA 1:100-1:2000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin 1:10-1:500 |

Images

Western Blot: PI 3-Kinase p85 alpha Antibody (S08-8H6) [NBP3-14962]
 - Western blot detection of PI 3-Kinase p85 alpha in HL-60, U2OS, C6 cell lysates using PI 3-Kinase p85 alpha Rabbit mAb (1:1000 diluted).
 Predicted band size: 84kDa. Observed band size: 84kDa.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-14962-100ul

| | |
|--------------------|--|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| H00005295-P01-10ug | Recombinant Human PI 3-Kinase p85 alpha GST (N-Term) Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-14962

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

