

Product Datasheet

Ferritin Heavy Chain Antibody (S03-3D9) NBP3-15028-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15028

Updated 11/17/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15028



NBP3-15028-100ul

Ferritin Heavy Chain Antibody (S03-3D9)

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S03-3D9
Preservative	0.01% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA
Target Molecular Weight	21 kDa


Product Description	
Host	Rabbit
Gene ID	2495
Gene Symbol	FTH1
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide of human Ferritin (Uniprot#: P02794)

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1/1000, Immunocytochemistry/ Immunofluorescence 1/20




Images


Western Blot: Ferritin Heavy Chain Antibody (S03-3D9) [NBP3-15028] - Western blot detection of Ferritin Heavy Chain in K562, CHO-K1, Hela cell lysates using Ferritin Heavy Chain Rabbit mAb (1:1000 diluted). Predicted band size: 21KDa. Observed band size: 21KDa.



Immunocytochemistry/Immunofluorescence: Ferritin Heavy Chain Antibody (S03-3D9) [NBP3-15028] - Immunofluorescence of Ferritin Heavy Chain (green) in Hela cells using Ferritin Rabbit mAb at dilution 1/5, and DAPI (blue)



Western Blot: Ferritin Heavy Chain Antibody (S03-3D9) [NBP3-15028] - Western blot detection of Ferritin in Raw264.7 cell lysates using Ferritin Heavy Chain Rabbit mAb (1:1000 diluted). Predicted band size: 21KDa. Observed band size: 21KDa.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15028-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBC1-18548	Recombinant Human Ferritin Heavy Chain Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15028

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

