

# Product Datasheet

## Sortilin Antibody (S02-6B4)

### NBP3-15095-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-15095](http://www.novusbio.com/NBP3-15095)

Updated 4/15/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-15095](http://www.novusbio.com/reviews/destination/NBP3-15095)



**NBP3-15095-100ul**


Sortilin Antibody (S02-6B4)

| Product Information         |                                                                                             |
|-----------------------------|---------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                      |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.      |
| Clonality                   | Monoclonal                                                                                  |
| Clone                       | S02-6B4                                                                                     |
| Preservative                | 0.01% Sodium Azide                                                                          |
| Isotype                     | IgG                                                                                         |
| Purity                      | Affinity purified                                                                           |
| Buffer                      | 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA                             |
| Target Molecular Weight     | 92 kDa                                                                                      |
| Product Description         |                                                                                             |
| Host                        | Rabbit                                                                                      |
| Gene Symbol                 | SORT1                                                                                       |
| Species                     | Human, Mouse, Rat                                                                           |
| Immunogen                   | Recombinant protein of human Sortilin (Uniprot#: Q99523)                                    |
| Product Application Details |                                                                                             |
| Applications                | Western Blot, ELISA, Immunohistochemistry-Paraffin                                          |
| Recommended Dilutions       | Western Blot, ELISA 1:100-1:2000, Immunohistochemistry-Paraffin                             |




## Images

Western Blot: Sortilin Antibody (S02-6B4) [NBP3-15095] - Western blot detection of Sortilin in Hela, A549, MCF-7 cell lysates using Sortilin Rabbit mAb (1:1000 diluted). Predicted band size: 92kDa. Observed band size: 92kDa.



Immunohistochemistry-Paraffin: Sortilin Antibody (S02-6B4) [NBP3-15095] - Immunohistochemistry of Sortilin in paraffin-embedded Human lung cancer tissue using Sortilin Rabbit mAb at dilution 1/50



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP3-15095-100ul**

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP1-89745PEP | Sortilin Recombinant Protein Antigen                |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-15095](http://www.novusbio.com/reviews/submit/NBP3-15095)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

