

Product Datasheet

Albumin Antibody (3O4L7)

NBP3-15415-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15415

Updated 7/31/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15415



NBP3-15415-100ul

Albumin Antibody (3O4L7)

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3O4L7
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol, 0.05% BSA
Target Molecular Weight	69 kDa

Product Description


Host	Rabbit
Gene ID	213
Gene Symbol	ALB
Species	Human, Mouse, Rat, Bovine
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of bovine Serum Albumin / BSA (P02769). KVASLRETYGDMADCCEKQEPERNECFLSHKDDSPDLPKLKPDPTLCDEFKA DEKKFWGKYLEIARRHPYFYAPELLYYANKYNGVVFQECQAEDKGAC

Product Application Details

Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:2000 - 1:10000, ELISA Recommended starting concentration is 1 ug/mL

Images

Western Blot: Albumin Antibody (3O4L7) [NBP3-15415] - Analysis of extracts of various cell lines, using Bovine Serum Albumin / BSA Rabbit mAb (NBP3-15415) at 1:10000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 1s.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15415-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-52084-0.05mg	Recombinant Human Albumin His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15415

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

