

Product Datasheet

Histone H1T Antibody - Azide and BSA Free NBP3-15574-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15574

Updated 10/2/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15574



NBP3-15574-100ul

Histone H1T Antibody - Azide and BSA Free

Product Information

Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Affinity purified
Buffer	PBS, 50% glycerol, pH7.3
Target Molecular Weight	22 kDa

Product Description

Host	Rabbit
Gene ID	3010
Gene Symbol	H1-6
Species	Mouse, Rat
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 110-209 of mouse Histone H1T (NP_034507.2). LSKKAASGNDKGKGKKSASAKAKKMGLPRASRSPKSSKTKAVKKPKATPTKAS GSGRKTGAKGVQQRKSPAKARAANPNSGKAKMVMQKTDLRKAAGRK


Product Application Details

Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500 - 1:2000, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:50 - 1:200




Images


Western Blot: Histone H1T Antibody [NBP3-15574] - Western blot analysis of extracts of various cell lines, using Histone H1T antibody (NBP3-15574) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry: Histone H1T Antibody [NBP3-15574] - Immunofluorescence analysis of mouse testis using Histone H1T Rabbit pAb (NBP3-15574) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunohistochemistry: Histone H1T Antibody [NBP3-15574] - Immunofluorescence analysis of rat testis using Histone H1T Rabbit pAb (NBP3-15574) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15574-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-56042PEP	Histone H1T Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15574

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

