

# Product Datasheet

## WDR57 Antibody - Azide and BSA Free NBP3-15611-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-15611](http://www.novusbio.com/NBP3-15611)

Updated 11/7/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-15611](http://www.novusbio.com/reviews/destination/NBP3-15611)



**NBP3-15611-100ul**

WDR57 Antibody - Azide and BSA Free

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 100 ul                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at -20C. Avoid freeze-thaw cycles.                                                    |
| Clonality           | Polyclonal                                                                                  |
| Preservative        | 0.01% Thimerosal                                                                            |
| Isotype             | IgG                                                                                         |
| Purity              | Affinity purified                                                                           |
| Buffer              | PBS, 50% glycerol, pH7.3                                                                    |


| Product Description |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Gene ID             | 9410                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Gene Symbol         | SNRNP40                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Species             | Human, Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Immunogen           | Recombinant fusion protein containing a sequence corresponding to amino acids 1-357 of human WDR57 (NP_004805.2).<br>MIEQQKRKGPELPLVPVKRQRHELLLGAGSGPGAGQQQATPGALLQAGPPRC<br>SSLQAPIMLLSGHEGEVYCCKFHPNGSTLASAGFDRILLWNVYGDCDNYATL<br>KGHSGAVMELHYNTDGSMLFSASTDKTVAVWDSETGERVKRLKGHTSFVN<br>SCYPARRGPQLVCTGSDDGTVKLWDIRKKAQITFQNTYQVLAVTFNDTSDQIISG<br>GIDNDIKVWDLRQNKLTYTMRGHADSVTGLSLSSSEGSYLLSNAMDNTVRVWD<br>VRPFAPKERCVKIFQGNVHNFENLLRCSWSPDGSKIAAGSADRFVYVWDTTS<br>RRILYKLPGHAGSINEVAFHPDEPIIIISASSDKRLYMGEIQ |

| Product Application Details |                             |
|-----------------------------|-----------------------------|
| Applications                | Western Blot                |
| Recommended Dilutions       | Western Blot 1:500 - 1:1000 |




## Images

Western Blot: WDR57 Antibody [NBP3-15611] - Western blot analysis of extracts of various cell lines, using WDR57 antibody (NBP3-15611) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Western Blot: WDR57 Antibody [NBP3-15611] - Western blot analysis of extracts of various cell lines, using WDR57 antibody (NBP3-15611) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP3-15611-100ul**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00009410-P01-10ug | Recombinant Human WDR57 GST (N-Term) Protein        |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-15611](http://www.novusbio.com/reviews/submit/NBP3-15611)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

