

Product Datasheet

Integrin alpha V/CD51 Antibody (10T9N2) NBP3-15647-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15647

Updated 7/31/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15647



NBP3-15647-100ul

Integrin alpha V/CD51 Antibody (10T9N2)


Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	10T9N2
Preservative	0.05% Proclin 300
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol, 0.05% BSA
Target Molecular Weight	116 kDa

Product Description	
Host	Rabbit
Gene ID	3685
Gene Symbol	ITGAV
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 949-1048 of human Integrin alpha V/CD51 (NP_002201.2). SLKSSASFNVIEFPYKNLPIEDITNSTLVTTNVTWGIQPAPMPVPVWVILAVLAG LLLLAVLVFVMYRMGFFKRVPPQEEQEREQLQPHENGEGNSET


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:1000 - 1:6000, ELISA Recommended starting concentration is 1 µg/mL, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/Immunofluorescence Validated from a verified customer review, Immunohistochemistry-Paraffin

Images


Western Blot: Integrin alpha V/CD51 Antibody (7O8N2) [NBP3-15647] - Analysis of extracts of various cell lines, using Integrin alpha V/CD51 (ITGAV/CD51) antibody (NBP3-15647) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.




ICC/IF: Integrin alpha V/CD51 Antibody (7O8N2) [NBP3-15647] -ITGAV Antibody [NBP3-15647] was used (1:100) for labeling and microscopy analysis of the distribution of ITGAV proteins in SUM159 cells. Validated customer review.




Immunohistochemistry: Integrin alpha V/CD51 Antibody (10T9N2) [Integrin alpha V/CD51] - Immunohistochemistry analysis of paraffin-embedded Mouse liver tissue using Integrin alpha V/CD51 Rabbit mAb at a dilution of 1:2000 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.



Immunohistochemistry: Integrin alpha V/CD51 Antibody (10T9N2) [Integrin alpha V/CD51] - Immunohistochemistry analysis of paraffin-embedded Human thyroid cancer tissue using Integrin alpha V/CD51 Rabbit mAb at a dilution of 1:2000 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.



Immunohistochemistry: Integrin alpha V/CD51 Antibody (10T9N2) [Integrin alpha V/CD51] - Immunohistochemistry analysis of paraffin-embedded Human spleen tissue using Integrin alpha V/CD51 Rabbit mAb at a dilution of 1:2000 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15647-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-85746PEP	Integrin alpha V/CD51 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15647

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

