

Product Datasheet

Smad4 Antibody (6B4D6)

NBP3-15673-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15673

Updated 1/29/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15673



NBP3-15673-100ul


Smad4 Antibody (6B4D6)

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	6B4D6
Preservative	0.05% Proclin 300
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol, 0.05% BSA
Target Molecular Weight	60 kDa
Product Description	
Host	Rabbit
Gene ID	4089
Gene Symbol	SMAD4
Species	Human, Mouse, Rat
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-552 of human Smad4 (NP_005350.1). PEYWCSIA YFEMDVQVGETFKVPSSCPIVTVDGYVDPSGGDRFCLGQLSNVH RTEAIERARLHIGKGVQLECKGEGDVWVRCLSDHAVFVQSYLDREAGRAPG DAVHKIYPSAYIKVFDLRQCHRQMQQQAATAQAAAAAQA AAVAGNIPGPGSVG GIAP AISLSAAAGIGVDDLRLRLCILRMSFVKGWGPDYPRQSIKETPCWIEIHLHR ALQLLDEV LHTMPIADPQPLD
Product Application Details	
Applications	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), Knockout Validated
Recommended Dilutions	Western Blot 1:1000 - 1:6000, ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements., Immunohistochemistry 1:500 - 1:2000, Immunoprecipitation 0.5ug-4ug antibody for 200µg-400µg extracts of whole cells, Immunohistochemistry-Paraffin 1:500 - 1:2000, Chromatin Immunoprecipitation (ChIP) 5ug antibody for 10µg-15µg of Chromatin, Knockout Validated




Images


Western Blot: Smad4 Antibody (1H7D3) [NBP3-15673] - Western blot analysis of extracts from normal (control) and Smad4 knockout (KO) 293T cells, using Smad4 antibody (NBP3-15673) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.




Western Blot: Smad4 Antibody (1H7D3) [NBP3-15673] -Analysis of various lysates using [KO Validated] Smad4 Rabbit mAb at 1:1000 dilution incubated overnight at 4c. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.




Western Blot: Smad4 Antibody (6B4D6) [Smad4] - Western blot analysis of various lysates using [KO Validated] Smad4 Rabbit mAb at 1:1000 dilution incubated overnight at 4C. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit . Exposure time: 3s.




Immunohistochemistry-Paraffin: Smad4 Antibody (1H7D3) [NBP3-15673] - Analysis of paraffin-embedded Human placenta tissue using [KO Validated] Smad4 Rabbit mAb at a dilution of 1:500 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.




Immunohistochemistry: Smad4 Antibody (6B4D6) [Smad4] - Immunohistochemistry analysis of paraffin-embedded Rat testis tissue using [KO Validated] Smad4 Rabbit mAb at a dilution of 1:500 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.




Immunohistochemistry: Smad4 Antibody (6B4D6) [Smad4] - Immunohistochemistry analysis of paraffin-embedded Mouse testis tissue using [KO Validated] Smad4 Rabbit mAb at a dilution of 1:500 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.




Chromatin Immunoprecipitation: Smad4 Antibody (6B4D6) [Smad4] - Chromatin immunoprecipitation analysis of extracts of HepG2 cells, using [KO Validated] Smad4 Rabbit mAb and rabbit IgG. The amount of immunoprecipitated DNA was checked by quantitative PCR. Histogram was constructed by the ratios of the immunoprecipitated DNA to the input.



Immunohistochemistry: Smad4 Antibody (6B4D6) [Smad4] - Immunohistochemistry analysis of paraffin-embedded Human cervix cancer tissue using [KO Validated] Smad4 Rabbit mAb at a dilution of 1:500 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.



Immunohistochemistry: Smad4 Antibody (6B4D6) [Smad4] - Immunohistochemistry analysis of paraffin-embedded Human placenta tissue using [KO Validated] Smad4 Rabbit mAb at a dilution of 1:500 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15673-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-24951PEP	Smad4 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15673

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

