

Product Datasheet

HuR/ELAVL1 Antibody (9Q6L3)

NBP3-15820-100ul

Unit Size: 100 ul

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-15820

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-15820



NBP3-15820-100ul

HuR/ELAVL1 Antibody (9Q6L3)

Product Information	
Unit Size	100 ul
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	9Q6L3
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.3), 50% glycerol, 0.05% BSA
Target Molecular Weight	36 kDa


Product Description	
Description	Novus Biologicals Knockout (KO) Validated Rabbit HuR/ELAVL1 Antibody (9Q6L3) (NBP3-15820) is a recombinant monoclonal antibody validated for use in IHC, WB, ELISA, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1994
Gene Symbol	ELAVL1
Species	Human, Mouse, Rat
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HuR/ELAVL1 (Q15717). MSGYEDHMAEDCRGDIGRTNLIVNYLPQNMTQDELRSLFSSIGEVESAKLIRD KVAGHSLGYGFVNYVTAKDAERAINTLNGLRLQSKTIKVSYPSS

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Knockout Validated
Recommended Dilutions	Western Blot 1:500 - 1:1000, ELISA, Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 1:50 - 1:200, Immunoprecipitation 0.5µg-4µg antibody for 200µg-400µg extracts of whole cells, Immunohistochemistry-Paraffin, Knockout Validated




Images


Western Blot: HuR/ELAVL1 Antibody (9Q6L3) [NBP3-15820] - Western blot analysis of extracts of various cell lines, using HuR/ELAVL1 antibody (NBP3-15820) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.




Western Blot: HuR/ELAVL1 Antibody (9Q6L3) [NBP3-15820] - Western blot analysis of extracts from normal (control) and HuR/ELAVL1 knockout (KO) 293T cells, using HuR/ELAVL1 antibody (NBP3-15820) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.




Immunohistochemistry: HuR/ELAVL1 Antibody (9Q6L3) [HuR/ELAVL1] - Immunohistochemistry analysis of paraffin-embedded Mouse spleen using [KO Validated] HuR/ELAVL1 Rabbit mAb at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Bufferr (pH 6.0) prior to IHC staining.




Immunohistochemistry: HuR/ELAVL1 Antibody (9Q6L3) [HuR/ELAVL1] - Immunohistochemistry analysis of paraffin-embedded Human kidney using [KO Validated] HuR/ELAVL1 Rabbit mAb at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Bufferr (pH 6.0) prior to IHC staining.




Immunocytochemistry/ Immunofluorescence: HuR/ELAVL1 Antibody (9Q6L3) [HuR/ELAVL1] - Confocal imaging of U-2 OS cells using [KO Validated] HuR/ELAVL1 Rabbit mAb . The cells were counterstained with alpha-Tubulin Mouse mAb (Green). DAPI was used for nuclear staining (blue). Objective: 100x.




Immunocytochemistry/ Immunofluorescence: HuR/ELAVL1 Antibody (9Q6L3) [HuR/ELAVL1] - Confocal imaging of NIH/3T3 cells using [KO Validated] HuR/ELAVL1 Rabbit mAb . The cells were counterstained with alpha-Tubulin Mouse mAb (Green). DAPI was used for nuclear staining (blue). Objective: 100x.




Immunohistochemistry: HuR/ELAVL1 Antibody (9Q6L3) [HuR/ELAVL1] - Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma using [KO Validated] HuR/ELAVL1 Rabbit mAb at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Bufferr (pH 6.0) prior to IHC staining.



Immunohistochemistry: HuR/ELAVL1 Antibody (9Q6L3) [HuR/ELAVL1] - Immunohistochemistry analysis of paraffin-embedded Human tonsil using [KO Validated] HuR/ELAVL1 Rabbit mAb at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Bufferr (pH 6.0) prior to IHC staining.



Immunoprecipitation: HuR/ELAVL1 Antibody (9Q6L3) [HuR/ELAVL1] - Immunoprecipitation analysis of 300 ug extracts of 293T cells using 3 ug HuR/ELAVL1 antibody . Western blot was performed from the immunoprecipitate using HuR/ELAVL1 antibody at a dilution of 1:1000.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-15820-100ul

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00001994-P01-2ug	Recombinant Human HuR/ELAVL1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-15820

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

