

# Product Datasheet

## GPRIN2 Antibody NBP3-17378-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-17378](http://www.novusbio.com/NBP3-17378)

Updated 5/3/2022 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-17378](http://www.novusbio.com/reviews/destination/NBP3-17378)



**NBP3-17378-100ul****GPRIN2 Antibody**

| <b>Product Information</b> |                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 100 ul                                                                                                           |
| <b>Concentration</b>       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>             | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| <b>Clonality</b>           | Polyclonal                                                                                                       |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                               |
| <b>Isotype</b>             | IgG                                                                                                              |
| <b>Purity</b>              | Affinity purified                                                                                                |
| <b>Buffer</b>              | PBS, pH 7.2, 40% glycerol                                                                                        |

| <b>Product Description</b> |                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>                | Rabbit                                                                                                                                                                         |
| <b>Gene ID</b>             | 9721                                                                                                                                                                           |
| <b>Gene Symbol</b>         | GPRIN2                                                                                                                                                                         |
| <b>Species</b>             | Human                                                                                                                                                                          |
| <b>Immunogen</b>           | This antibody was developed against Recombinant Protein corresponding to amino acids:<br>NVSTMGGGDLCLRLAPSAAAMQRSHSDLVRSTQMRGHSGARKASLSCSALG<br>SSPVHRAQLQPGGTSGQGGQAPAGLERDLA |

| <b>Product Application Details</b> |                                                                                                        |
|------------------------------------|--------------------------------------------------------------------------------------------------------|
| <b>Applications</b>                | Immunocytochemistry/Immunofluorescence                                                                 |
| <b>Recommended Dilutions</b>       | Immunocytochemistry/Immunofluorescence 0.25-2 ug/ml                                                    |
| <b>Application Notes</b>           | Immunocytochemistry/Immunofluorescence, PFA/Triton X-100 is recommended for fixation/permeabilization. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-17378](http://www.novusbio.com/reviews/submit/NBP3-17378)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

