

# Product Datasheet

## ODF3L2 Antibody NBP3-17406-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-17406](http://www.novusbio.com/NBP3-17406)

Updated 5/3/2022 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-17406](http://www.novusbio.com/reviews/destination/NBP3-17406)



**NBP3-17406-100ul**

ODF3L2 Antibody

| Product Information         |                                                                                                                                                                        |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                                                 |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                       |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                                                                 |
| Clonality                   | Polyclonal                                                                                                                                                             |
| Preservative                | 0.02% Sodium Azide                                                                                                                                                     |
| Isotype                     | IgG                                                                                                                                                                    |
| Purity                      | Affinity purified                                                                                                                                                      |
| Buffer                      | PBS, pH 7.2, 40% glycerol                                                                                                                                              |
| Product Description         |                                                                                                                                                                        |
| Host                        | Rabbit                                                                                                                                                                 |
| Gene ID                     | 284451                                                                                                                                                                 |
| Gene Symbol                 | ODF3L2                                                                                                                                                                 |
| Species                     | Human                                                                                                                                                                  |
| Immunogen                   | This antibody was developed against Recombinant Protein corresponding to amino acids:<br>EIPGPGQYDSPDANTYRQRLPAFTMLGRPRAPRPLEETPGPGAHCPEQVTVN<br>KARAPAFSMGIRHASKRASTM |
| Product Application Details |                                                                                                                                                                        |
| Applications                | Immunocytochemistry/Immunofluorescence                                                                                                                                 |
| Recommended Dilutions       | Immunocytochemistry/Immunofluorescence 0.25-2 ug/ml                                                                                                                    |
| Application Notes           | Immunocytochemistry/Immunofluorescence, PFA/Triton X-100 is recommended for fixation/permeabilization.                                                                 |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-17406](http://www.novusbio.com/reviews/submit/NBP3-17406)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

