

# Product Datasheet

## IL1RAPL1 Antibody (B-L43) [PE] NBP3-18110PE

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-18110PE](http://www.novusbio.com/NBP3-18110PE)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-18110PE](http://www.novusbio.com/reviews/destination/NBP3-18110PE)



**NBP3-18110PE**

IL1RAPL1 Antibody (B-L43) [PE]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | B-L43                                                                                       |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1 Kappa                                                                                  |
| Conjugate           | PE                                                                                          |
| Purity              | Affinity purified                                                                           |
| Buffer              | PBS                                                                                         |

| Product Description     |                                                                                  |
|-------------------------|----------------------------------------------------------------------------------|
| Host                    | Mouse                                                                            |
| Gene ID                 | 11141                                                                            |
| Gene Symbol             | IL1RAPL1                                                                         |
| Species                 | Human                                                                            |
| Specificity/Sensitivity | Recognizes natural and recombinant human IL1RAPL1                                |
| Immunogen               | Recombinant human IL1RAPL1. Hybridoma: Myeloma X63/AG.8653 x Balb/c spleen cells |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Flow Cytometry                                                         |
| Recommended Dilutions       | Flow Cytometry                                                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-18110PE**

---

|                    |                                                    |
|--------------------|----------------------------------------------------|
| NBP1-43319PE-0.5ml | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [PE] |
| 201-LB-005         | IL-1 beta/IL-1F2 [Unconjugated]                    |
| 9949-MR-050        | IL1RAPL1 [Unconjugated]                            |
| NB100-56566        | TLR4 Antibody (76B357.1) - BSA Free                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-18110PE](http://www.novusbio.com/reviews/submit/NBP3-18110PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

