

# Product Datasheet

## PEX10 Antibody Blocking Peptide NBP3-18137PEP

Unit Size: 0.05 mg

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-18137PEP](http://www.novusbio.com/NBP3-18137PEP)

Updated 6/30/2022 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-18137PEP](http://www.novusbio.com/reviews/destination/NBP3-18137PEP)



**NBP3-18137PEP****PEX10 Antibody Blocking Peptide**

| <b>Product Information</b>         |                                                                                                                                                             |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.05 mg                                                                                                                                                     |
| <b>Concentration</b>               | 0.2 mg/ml                                                                                                                                                   |
| <b>Storage</b>                     | Store at -20C. Avoid freeze-thaw cycles.                                                                                                                    |
| <b>Preservative</b>                | 0.02% Sodium Azide                                                                                                                                          |
| <b>Purity</b>                      | N/A                                                                                                                                                         |
| <b>Buffer</b>                      | PBS pH 7.2 (10 mM NaH <sub>2</sub> PO <sub>4</sub> , 10 mM Na <sub>2</sub> HPO <sub>4</sub> , 130 mM NaCl) containing 0.1% bovine serum albumin             |
| <b>Product Description</b>         |                                                                                                                                                             |
| <b>Description</b>                 | <p>A blocking peptide containing 10 amino acids near the carboxy terminus of PEX10.</p> <p><b>Source:</b> <i>Synthetic</i></p> <p>(Accession #: O60683)</p> |
| <b>Gene ID</b>                     | 5192                                                                                                                                                        |
| <b>Gene Symbol</b>                 | PEX10                                                                                                                                                       |
| <b>Species</b>                     | Human                                                                                                                                                       |
| <b>Product Application Details</b> |                                                                                                                                                             |
| <b>Applications</b>                | Antibody Competition                                                                                                                                        |
| <b>Recommended Dilutions</b>       | Antibody Competition                                                                                                                                        |
| <b>Application Notes</b>           | This peptide is useful as a blocking peptide for NBP3-18137. For further blocking peptide related protocol, click <a href="#">here</a> .                    |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NBP3-18137PEP**

---

|                    |                                                  |
|--------------------|--------------------------------------------------|
| H00005192-P01-10ug | Recombinant Human PEX10 GST (N-Term) Protein     |
| MAB1455            | Albumin Antibody (188835) [Unconjugated] - Serum |
| NBP3-18137         | PEX10 Antibody - BSA Free                        |
| 2914-HT-100MG      | Holo-Transferrin                                 |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-18137PEP](http://www.novusbio.com/reviews/submit/NBP3-18137PEP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

