

Product Datasheet

PEX3 Antibody - BSA Free NBP3-18138

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-18138

Updated 6/13/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-18138



NBP3-18138

PEX3 Antibody - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS

Product Description	
Host	Rabbit
Gene ID	8504
Gene Symbol	PEX3
Species	Human, Mouse, Rat
Immunogen	Antibody was raised against a peptide corresponding to 12 amino acids near the carboxyl terminal of human PEX3. The immunogen is located within last 50 amino acids of PEX3.


Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 2 ug/mL, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence 20 ug/mL

Images


Western Blot: PEX3 Antibody - BSA Free [NBP3-18138] - Western Blot Validation in Rat Tissues. Loading: 15 ug of lysates per lane. Antibodies: PEX3 NBP3-18138, 2 ug/mL, 1h incubation at RT in 5% NFDM/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10,000 dilution.




Immunocytochemistry/Immunofluorescence: PEX3 Antibody - BSA Free [NBP3-18138] - Immunofluorescent analysis of 4% paraformaldehyde-fixed HeLa cells labeling PEX3 with NBP3-18138 at 20 ug/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (green) and DAPI staining (blue).




Immunohistochemistry: PEX3 Antibody - BSA Free [NBP3-18138] - Immunofluorescent analysis of 4% paraformaldehyde-fixed mouse thymus labeling PEX3 with NBP3-18138 at 20 ug/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (green) and DAPI staining (blue).




Western Blot: PEX3 Antibody - BSA Free [NBP3-18138] - WB Validation in Human Cell Lines. Loading: 15 ug of lysate Antibodies: PEX3 NBP3-18138, 2 ug/mL, 1 h incubation at RT in 5% NFD/MTBST. Secondary: Goat Anti-Rabbit IgG HRP conjugate at 1:10,000 dilution.




Western Blot: PEX3 Antibody - BSA Free [NBP3-18138] - Western Blot Validation in Mouse Tissues. Loading: 15 ug of lysates per lane. Antibodies: PEX3 NBP3-18138, 2 ug/mL, 1h incubation at RT in 5% NFD/MTBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10,000 dilution.



Immunohistochemistry: PEX3 Antibody - BSA Free [NBP3-18138] - Immunofluorescent analysis of 4% paraformaldehyde-fixed human spleen labeling PEX3 with NBP3-18138 at 20 ug/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (green) and DAPI staining (blue).



Immunohistochemistry: PEX3 Antibody - BSA Free [NBP3-18138] - Immunofluorescent analysis of 4% paraformaldehyde-fixed rat testis labeling PEX3 with NBP3-18138 at 20 ug/mL, followed by goat anti-rabbit IgG secondary antibody at 1/500 dilution (green) and DAPI staining (blue).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP3-18138

NBP3-18138PEP	PEX3 Antibody Blocking Peptide
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7156	Goat anti-Rabbit IgG (H+L) Secondary Antibody
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-18138

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

