

Product Datasheet

YIF1B Antibody - BSA Free

NBP3-18392

Unit Size: 100 ug

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-18392

Updated 3/18/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-18392



NBP3-18392

YIF1B Antibody - BSA Free

Product Information

| | |
|----------------------|---------------------------------|
| Unit Size | 100 ug |
| Concentration | 1 mg/ml |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Affinity purified |
| Buffer | Tris-citrate/phosphate pH 7 -8. |

Product Description


| | |
|--------------------|--|
| Host | Rabbit |
| Gene ID | 90522 |
| Gene Symbol | YIF1B |
| Species | Human |
| Immunogen | Recognizes region between residues Between between 264 and 314.. |

Product Application Details

| | |
|------------------------------|--|
| Applications | Immunoprecipitation |
| Recommended Dilutions | Immunoprecipitation 15-20 ug/mg lysate |

Images

Immunoprecipitation:[NBP3-18392]-Detection of human YIF1B by western blot of immunoprecipitates. Samples: Whole cell lysate from K-562 cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-YIF1B antibody used for IP at 6 µg per reaction. YIF1B was also immunoprecipitated by rabbit anti-YIF1B antibody. For blotting immunoprecipitated YIF1B, was used at 0.4 mg/ml. Detection: Chemiluminescence with an exposure time of 3 seconds.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-18392

| | |
|------------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP2-51700-0.1mg | Recombinant Human YIF1B His Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-18392

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

