

Product Datasheet

CD9 Antibody (HI9a) - BSA Free

NBP3-18795

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP3-18795

Updated 3/3/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP3-18795



NBP3-18795

CD9 Antibody (HI9a) - BSA Free

Product Information

Unit Size	100 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	HI9a
Preservative	0.09% Sodium Azide
Isotype	IgG1 Kappa
Purity	Affinity purified
Buffer	PBS, pH 7.2

Product Description


Description	May be stored at 4C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20C. Aliquots are stable for 12 months after receipt. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.
Host	Mouse
Gene ID	928
Gene Symbol	CD9
Species	Human, Bovine, Canine, Equine, Monkey, Rabbit, Sheep
Immunogen	A recombinant protein of human CD9
Notes	May be stored at 4C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20C. Aliquots are stable for 12 months after receipt. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

Product Application Details

Applications	Flow Cytometry
Recommended Dilutions	Flow Cytometry 0.5 ug/10 million cells in 100 ul
Application Notes	For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

Images

Flow Cytometry: CD9 Antibody (HI9a) [NBP3-18795] - Human platelets stained with NBP3-18795, then detected with anti-mouse IgGs FITC.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP3-18795

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
236-EG-200	EGF [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP3-18795

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

