

# Product Datasheet

## RAB8A Antibody (S06-8A7)

### NBP3-19493-100ul

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP3-19493](http://www.novusbio.com/NBP3-19493)

Updated 4/15/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP3-19493](http://www.novusbio.com/reviews/destination/NBP3-19493)



**NBP3-19493-100ul**


RAB8A Antibody (S06-8A7)

| Product Information         |                                                                                                                            |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                     |
| Concentration               | 0.3 mg/ml                                                                                                                  |
| Storage                     | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                                     |
| Clonality                   | Monoclonal                                                                                                                 |
| Clone                       | S06-8A7                                                                                                                    |
| Preservative                | 0.01% Sodium Azide                                                                                                         |
| Isotype                     | IgG                                                                                                                        |
| Purity                      | Affinity purified                                                                                                          |
| Buffer                      | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.05% BSA                                                             |
| Target Molecular Weight     | 24 kDa                                                                                                                     |
| Product Description         |                                                                                                                            |
| Host                        | Rabbit                                                                                                                     |
| Gene Symbol                 | RAB8A                                                                                                                      |
| Species                     | Human, Mouse, Rat                                                                                                          |
| Immunogen                   | A synthetic peptide of human RAB8A (Uniprot # P61006)                                                                      |
| Product Application Details |                                                                                                                            |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                 |
| Recommended Dilutions       | Western Blot, Immunohistochemistry 1:50-1:200, Immunocytochemistry/ Immunofluorescence 1:50, Immunohistochemistry-Paraffin |




## Images


**Western Blot: RAB8A Antibody (S06-8A7) [NBP3-19493]** - Western blot detection of RAB8A in K562, C6, 3T3, HeLa cell lysates using NBP3-19493 (1:1000 diluted). Predicted band size: 24kDa. Observed band size: 24kDa.



**Immunocytochemistry/Immunofluorescence: RAB8A Antibody (S06-8A7) [NBP3-19493]** - Immunofluorescence of RAB8A (green) in HeLa cells using NBP3-19493 at dilution 1/200, and DAPI (blue)



**Immunohistochemistry-Paraffin: RAB8A Antibody (S06-8A7) [NBP3-19493]** - Immunohistochemistry of RAB8A in paraffin-embedded Human colon cancer tissue using NBP3-19493 at dilution 1/1





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP3-19493-100ul**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00004218-P01-10ug | Recombinant Human RAB8A GST (N-Term) Protein        |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP3-19493](http://www.novusbio.com/reviews/submit/NBP3-19493)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

